



medava[®]
SOFTTOUCH GLOVES

NITRILE EXAMINATION GLOVES

PRODUCT SPECIFICATION

- CHEMO
- NON STERILE
- POWDER FREE
- LATEX FREE

Grade	Chemo Grade	Surface	Smooth	Glove Length(mm)	≥230
Material	NBR	Packing	200pcs/box,10boxes/ctn	Weight(g)	3.5±0.3g(Medium Size)
Cuff	Beaded	Shelf-life	5 years	Powder Level	≤2mg/glove
		SIZE			
		Small	Medium	Large	Extra Large
Weight(g)		3.2±0.3g	3.5±0.3g	3.8±0.3g	4.1±0.3g
Physical Dimension	Length (mm)	≥230	≥230	≥230	≥230
	Palm Width (mm)	80±10	95±10	110±10	120±10
	Cuff Thickness(mm)	0.07±0.03	0.07±0.03	0.07±0.03	0.07±0.03
	PalM Thickness(mm)	0.08±0.03	0.08±0.03	0.08±0.03	0.08±0.03
	Finger Thickness(mm)	0.10±0.03	0.10±0.03	0.10±0.03	0.10±0.03
Physical Property	Tensile Strength(MPa)	Before Aging	16		
		After Aging	14		
	Elongation(%)	Before Aging	500		
		Before Aging	400		
Water Test Before AQL		AQL1.5	AQL1.5	AQL1.5	AQL1.5
Appearance AQL		Critical AQL1.0 Major AQL2.5 Minor 4.0	Critical AQL1.0 Major AQL2.5 Minor 4.0	Critical AQL1.0 Major AQL2.5 Minor 4.0	Critical AQL1.0 Major AQL2.5 Minor 4.0
Certifications: 510K , ISO9001 ISO13485					
Sampling Procedures according to the standard of ISO13485, ASTM D6319, ASTM D6978, ASTM F1671 Inspection Level for physical dimensions and force /elongation at break is according to S-2 . water test according to G-1					

	XS	Small	Medium	Large	Extra large
CASE SKU	AMCC2000-XS	AMCC2000-S	AMCC2000-M	AMCC2000-L	AMCC2000-XL

MEDAVA® Authorized Distributor



1487 Trae Lane
Lithia Springs, GA 30122
www.unitemedical.net

P/F: 833.TO-Unite
(833.868.6483)
Email: Orders@UniteMedical.net