



ANR Technologies Pte Ltd

Chamfering Pliers

SPECIFICATION

The manual disc cutter is used for coin cell electrode with customized cutting die. Simple operation and light in weight. The cut-out discs can be used as electrode/seperator for battery research.



DETAILS

Category	Value
SKU	CNA-7-9 for R3 R5 R10
Shape	Circular
Diameter	R3,R5,R10
Cutting Thickness	Maximum thickness: 0.3mm for soft metallic foil
Packaging Dimensions	18*9*2.5cm
Weight	0.4kg

