



DNTLworks Equipment Corporation

Specifications

2400.02.00

ProCart III Self-Contained Treatment Console 120V



Amperage	10 Amps. 12 Amps Max with Accessories
Compressor Type	Compressor and Vacuum Pump: Oil Free
Dimensions	24"x26"x32"
Element Filtration	50 micron Vacuum System
Element Filtration	5 micron Air System
Flow Rate	5.0scfm @ 10inhg. Vacuum
Flow Rate	2.5scfm @ 40psi. Compressor
Font	Adobe Garamond Pro Regulator
Hertz	60hz
Horse Power	1/3Hp Compressor/Vacuum
Mic.	Refining Cartridge About \$90.00
Operating Pressure	70 to 100 psi Compressor
Regulated Pressure	40psi Water System
Regulated Pressure	65psi Air System
Shipping Size	30"x27"x40"
Shipping Weight	140Lbs
Sound Level	44db @3'
Voltage	120Vac
Volume Air Tank	1 gal, 3.7 Liter
Volume Vacuum Container	2qt 1.9Liter
Volume Water Container	2 Each, 1 Liter Bottle
Wattage	1200 to 1440
Weight	113 LBS