

# **Mobile Evacuation Unit**

The PC 2530 Mobile Vac System is perfect where central installation is impractical.

# Moves with Ease

Reinforced large polyethylene carpet casters Efficient

Variable controlled motor speed can evacuate at the rate of 75 cfm, or my be adjusted to the user's preference

# Simple Maintenance

Non-corrosive collection canister cleans easily and is equipped with an automatic overflow shut-off

# Perfect Backup

Invaluable backup when the central system breaks down. The user just rolls it out, plugs it in, and continues operating.

# **Cost Effective**

Low cost leader of oral evacuation units. Excellent for limited budgets.

# Compact

Measures only 14" x 12" x 31". Easily placed at chair side and stores inconspicuously.

The PC 2530 Mobile Vac System is a portable, self contained vacuum system suitable for use by the dentist in a temporary location or when it is not feasible to install a central vacuum system.

# Package include the following:

- Equipped with a muffler system for quiet operation
- Hve and SE valves on each suction hose
- Activates only when Hve or SE valves are lifted from hanger
- High volume suction
- White high impact resistant polyethylene case
- 1 year warranty



Mobile Vacuum System



# **Specifications**

Electrical 7/8 HP 115 VAC, 50/60 Hz 230 VAC, 50/60 Hz 4.0 Amps



Packing Dimensions (L x W x H) 17" x 16" x 33' 70 lbs

TPC Advanced Technology, City of Industry, California • Tel: (626) 810-4337 • Fax: (626) 810-4245 tpcdental.com