



DNTLworks Equipment Corporation

Specification

2600.34.00

ProCart I Self-Contained Treatment Cart 120V



Amperage	10 Amps. 12 Amps Max with Accessories
Compressor Type	Compressor and Vacuum Pump: Oil-Free
Dimensions	23"x23"x32 to 42"
Flow Rate	2.5scfm @ 40psi Compressor
Flow Rate	4.8scfm @ 10 Inhg Vacuum
Hertz	60hz
Horse Power	1/3Hp Compressor/Vacuum
Operating Pressure	70 to 100psi Compressor:
Regulated Pressure	65psi Air Sysytem
Regulated Pressure	40psi Water System:
Shipping Size	27x27x20" Both Delivery, Compressor
Shipping Weight	59lbs Delivery, 74lbs Compressor
Sound Level	51db @ 3'
Voltage	120Vac
Volume Air Tank	3/4gal
Volume Vacuum Container	2qt, 1.9 Liter
Volume Water Container	1 Liter
Wattage	1200 to 1440
Weight	96 Lbs