




CLEANAIR PRO

“Filtration for your vacuum system offers a cleaner environment from viruses and bacteria.”

- ✓ ULPA filters are recommended in applications where extremely hazardous particles can pose a threat to the operator and patient
- ✓ ULPA filters produce 99.995% particle removal at 0.12 micron
- ✓ ULPA filters are superior to a HEPA filter which is 99.97% at 0.3 micron
- ✓ Our ULPA filters are used in clean rooms
- ✓ Can be adapted to any DNTLworks equipment without reduction of vacuum efficiency
- ✓ Filtration system comes with everything you need to adapt to DNTLworks equipment
- ✓ Why take the chance when it comes to you and your patient's health
- ✓ Buy an American made DNTLworks filtration system today!



Materials	Specifications
Element Filtration	99.9995% @ 0.12 Micron
Filter Size	14x7.5x9
Filter Weight	5.0lbs
Filter Material	ULPA
Vacuum Flow	Max: 8scfm