



DNTLworks Equipment Corporation

Specifications

4017.03.00

Basic Aluminum Patient Chair with Hydraulic Base



Backrest Adjustment	90 to 180 degrees supine
Dimensions	75x22-1/2 Adjustable Height 19" to 31"
Folded Dimensions	60-1/2"x 22-1/2"x12"
Maximum Load Weight	500 Lbs
Shipping Size	55"x25"x16" Chair, 27"x27"x20" Base
Shipping Weight	54 Lbs Chair, 56 Lbs Base
Weight	76 Lbs Total
Weight	34 Lbs chair
Weight	42 Lbs hydraulic base