



DNTLworks Equipment Corporation

Specifications

4025.03.00

Supreme Aluminum Patient Chair with Hydraulic Base



Backrest Adjustment	90 to 180 degrees supine
Dimensions	74x22-1/2, 18 to 30" Seat Height.
Folded Dimensions	58"x52"
Maximum Load Weight	500 Lbs
Shipping Size	55"x25"x16" Chair, 27"x27"x20" Base
Shipping Weight	56 Lbs Chair, 57 Lbs Base
Weight	77 lbs, 500lbs Max
Weight	35 Lbs Chair
Weight	77 Lbs Total