



# DNTLworks Equipment Corporation

## Specifications

6238.00.00

ProBrite Portable LED 10 Light W/Swing Arm



Amperage	1mA
Focal Area	At 24"/Foot Candles: 6-1/4"/1210fc
Focal Area	At 18"/Foot Candles: 4-5/8"/2180fc
Focal Area	At 12"/Foot Candles: 3-1/4"/4882fc
Goose Neck	28"
Lux	6,500 to 45,000
Shipping Size	33"x18"x10"
Shipping Weight	15 Lbs
Voltage	100 to 240Vac
Weight	8 Lbs