



# DNTLworks Equipment Corporation

## Specifications

6336.00.00

ProBrite Portable LED 3 Light W/Surface Mount



Amperage	.7mA
Focal Area	At 12"/Foot Candles:4"/4890fc
Focal Area	At 18"/Foot Candles: 5"/2110fc
Focal Area	At 24"/Foot Candles: 6"/1240fc
Goose Neck	28"
Lux	42,000
Shipping Size	33"x18"x10"
Shipping Weight	12 Lbs
Voltage	100 to 240Vac
Weight	5 Lbs