

# FORZA™

Clinical Electric Systems



# FORZA ELS™

*Muscular yet Mindful...*

The Forza ELS sets the standard for clinical micromotors with a powerful 3.0 Ncm of torque in a compact, lightweight design. Smooth, quiet, and versatile, the Forza ELS features two programmable settings and an easy-to-read digital speed display. Whether cutting through obsolete crowns or refining margins, the Forza ELS delivers the power and control demanded by dentists today.

## Forza ELS Electric Micromotor

Item Number: 5021306U0

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.0 Ncm of Torque

76 Grams Motor Weight

3 Year Warranty

