

FORZA™

Clinical Electric Systems



FORZA ELM™

Changing the Way You Think About Electric Handpieces...

The Forza ELM is arguably the most advanced clinical micromotor. The Forza ELM is a celebration of details, from the lightweight titanium motor to the LED fiber optics. With the lightest motor in its class and a powerful 3.4 Ncm of torque, this micromotor can be integrated into any dental delivery system. The Forza ELM will change the way you think about electric handpieces.

Forza ELM Electric Micromotor

Item Number: 5021312U0

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.4 Ncm of Torque

61 Grams Motor Weight

3 Year Warranty

LED Optics



FORZA™

Clinical Electric Systems



FORZA ELS™

Muscular yet Mindful...

The Forza ELS sets the standard for clinical micromotors with a powerful 3.0 Ncm of torque in a compact, lightweight design. Smooth, quiet, and versatile, the Forza ELS features two programmable settings and an easy-to-read digital speed display. Whether cutting through obsolete crowns or refining margins, the Forza ELS delivers the power and control demanded by dentists today.

Forza ELS Electric Micromotor

Item Number:5021306U0

System Includes Motor, Cord, Control Unit and AC Adapter

Two Programmable Settings

Motor Speed Range: 2,000-40,000 rpm

Attachments Provide a System Speed Range of 100-200,000 rpm

3.0 Ncm of Torque

76 Grams Motor Weight

3 Year Warranty

