


GPRO® F1

Gpro F1 adopts wireless, vertical type, and single button design, bringing easier obturation operation.



Quick Heating within 15 Seconds

Reduces time for waiting and increase your working efficiency.



Safer Protection Mechanism

Automatic power-off after working for 10min.



Ergonomic Design

Easy to hold and operate. Unique double displays design enables convenient operations in different angles and teeth positions.



360° Rotatable Injection Needle

Injection needle made by special technology can be rotated freely to facilitate multi-angle filling. Easy to reach the apical 3mm of canal brings better filling.

Key Specifications

Size	31.9mmx152.5mmx114.9mm
Weight	170g
Charging Base Weight	207g
Charging Base Size	75.5mmx149.7mmx62.6mm
Classification	Class II(AC/DC power adapter)
Optional preset temperatures	150°-180°-200°-230°
Time consumption for charging	About 2.5h
Power Input	AC100V-240V 50/60Hz 0.8mA
Power Output	DC15V/1.6A
Battery capacity(Chargeable battery)	2000mAh
Warranty	12 Months

Order Information



Gutta-percha Obturation Filling Machine

Code: GP2002
Model: GPRO F1



Gutta-percha Obturation System

Code: GP2001
Model: GPRO S1



Beyes Dental Canada Inc.
23-595 Middlefield Road
Toronto, Ontario, M1B 3M5
Canada
Tel: 1-855-603-1888
Email: info@beyes.ca
Web: www.beyes.ca



Lotus NL B.V.
Koningin Julianaplein 10, le Verd,
2595AA, The Hague,
Netherlands
Tel: +31645171879 (English)
+31626669008 (Dutch)
Email: info@beyes.ca

Your Dealer for Beyes Products

Scan the QR code
for more products



Printed in Canada
File #: ENB034
Rev.1/22.01.21

© Beyes Canada Inc 2021. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance
with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under Canadian
Medical Devices Conformity Assessment System (CMDCAS).