



Model #L690-LED  
Ceiling Mounted Lucent LED  
Operatory Light



### Features:

- light head rotation
- The contoured handles are conveniently located for operator positioning and asepsis protection
- 3 step dimmer mode
- Post mounted, Ceiling Mounted and Wall mount options
- Removable plastic splash shield
- Optional: motion sensor
- 2 year warranty
- Packing dimensions L600-LED: 38 x 17 x 15 28 lbs.

### Specifications:

Absorbed Power:	Max 8vA
Input voltage:	110V/220V, 60/50 Hz
Output voltage:	24 Volts AC
Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Cooling System:	Static Heat sink (no fan)
Illumination area:	100mm x 180mm @ 700mm

