

	No	Type		No	Type		No	Type
<b>AB Dental</b>			<b>BEGO</b>			<b>Brain Base</b>		
Narrow Platform (NP)	55021	11	Semados Mini	55012	2	Mytis Arrow B-Type 3.7-4.6	55014	4
Internal Hex Platform (SP)	55049	41	Semados S / RI 3.25-5.5	55034	26	Mytis Arrow BS-type 3.7-4.6	55014	4
Conical Platform (CP)	55046	38	Semados RS / RSX 3.0	55082	74	Mytis Arrow C-Type 3.3	55014	4
			Semados RS / RSX 3.75-5.5	55034	26	Mytis Arrow C-Type 3.7-4.6	55020	10
<b>ADIN</b>			Semados SC / SCX 3.25-5.5	55034	26	Mytis Arrow E-Type 3.3-4.6	55014	4
Touareg S/OS	55034	26				Mytis Arrow Me-Type 4.0-5.0	55020	10
Touareg Closefit UNP 2.75	55059	51	<b>Bio3</b>			Mytis Arrow St-Type 4.0-5.0	55020	10
Touareg Closefit NP 3.0	55059	51	Advanced and Progressive 3.3-3.8	55025	16			
Touareg Closefit RP 3.5	55046	38	Advanced and Progressive 4.2-5.0	55067	59	<b>Bredent</b>		
Touareg CloseFit WP 4.3, 5.0	55020	10				SKY	55025	16
Swell 3.3	55034	26	<b>BioHorizons</b>			copaSKY	55027	19
			Laser-Lok 3.0 Grey	55034	26	mini2SKY	55033	25
<b>AIRES</b>			Internal Laser-Lok 3.5 Yellow	55049	41	fast & fixed Abutment	55083	75
3.3	55049	41	Internal Laser-Lok 4.5 Green	55051	43			
			Internal Laser-Lok 5.7 Blue	55051	43	<b>BTI</b>		
<b>Alfa Gate</b>			External 4.0 Platform	55011	1	Interna Universal	55051	43
Bioactive, All diam	55034	26				Interna Wide	55051	43
MAX, All diam	55034	26	<b>BioMate</b>					
M+ NP 3.5	55044	36	BioMate 3.3	55074	66	<b>BTK</b>		
M+ RP 4.3 - 5.0	55037	29	BioMate 4.8	55081	73	External Hex EN	55038	30
Porous, All diam	55034	26				External Hex ER	55011	1
S-Line 3.0	55059	51	<b>Biomet/3i</b>			External Hex EW	55011	1
			External Hex	55011	1	Internal Hex IR	55049	41
<b>Alpha Dent</b>			Certain internal Purple	55036	28	Internal Hex IM	55051	43
Conical connection Narrow Platform	55025	16	Certain internal Purple Platform 3.4	55036	28	Internal Hex IW	55051	43
Conical connection Wide Platform	55057	49	Certain internal Blue, Yellow, Green	55066	58	Hexagonal Morse Taper KR	55081	73
Hexagon connection Narrow Platform	55039	31	T3 Purple	55036	28	Hexagonal Morse Taper KW	55056	48
Hexagon connection Regular Platform	55057	49	T3 Blue, Yellow, Green	55066	58	Octagonal Morse Taper SR	55014	4
<b>Alpha-Bio</b>			<b>Bioner</b>			<b>C-Tech</b>		
Internal Hex	55034	26	HIKELT 3.8-4.7	55034	26	EL	55016	6
Conical Hex 3.2	55065	57	Top DM 3.5-5.0	55046	38			
Conical Hex 3.75-5.0	55051	43				<b>Camlog</b>		
			<b>Bionika</b>			Screw-Line, Root-Line, Progressive 3.3	55036	28
<b>Anker</b>			Cortical 4.0	55034	26	Screw-Line, Root-Line, Progressive 3.8, 4.3	55031	23
SB III 3.5	55023	14	CortiLog (CCL) 3.3	55065	57	Screw-Line, Root-Line, Progressive 5.0, 6.0	55032	24
SB III 4.0-5.0	55020	10	CortiLog (CCL) 3.8	55023	14	iSy	55065	57
			CortiLog (CCL) 5.0	55075	67			
<b>Ankylos</b>			<b>Bionnovation</b>			<b>Ceralog</b>		
C/X	55025	16	Biomorse	55023	14	Hexalobe ceramic	55076	68
<b>Anthogyr</b>			<b>Biotech Dental/Ids</b>			<b>CLC Scientific</b>		
Axiom Bone Level 3.4 - 5.2	55065	57	Kontakt, Kontakt-S 3.6 - 5.4	55045	37	Conic 3,5-6,0	55081	73
Axiom Tri-lobe	55065	57	Kontakt-S+ 4.5-5.5	55045	37			
Axiom Tissue Level Narrow, Regular	55050	42	Reflect Aspire	55081	73	<b>Conelog</b>		
			Reflect Certus	55019	9	Screw-Line, Progressive 3.3	55066	58
<b>Argon</b>			Reflect Tapered	55034	26	Screw-Line, Progressive 3.8, 4.3	55066	58
K3pro Rapid 3.0 - 3.5	55026	17	Reflect Recover	55020	10	Screw-Line, Progressive 5.0	55075	67
K3pro Rapid Platf 3mm	55019	9	Reflect Rapid	55020	10			
K3pro Rapid Platf 3mm Short	55044	36				<b>Cowellmedi</b>		
			<b>Biotem</b>			Inno Int Octa	55014	4
<b>Astra Tech</b>			AR 3.7 (Yellow M)	55066	58	Inno Sub Hex	55037	29
Osseospeed EV 3.0	55058	50	AR 4.2-4.6 (Green R)	55037	29	<b>Dentack</b>		
Osseospeed EV 3.6	55046	38	AR 5.1 (Blue W)	55037	29	Internal Hex	55056	48
Osseospeed EV 4.2	55057	49	BR 3.3 (Yellow M)	55012	2			
Osseospeed EV 4.8	55055	47	BR 3.75-4.0 (Green R)	55014	4	<b>Dental Tech</b>		
Osseospeed EV 5.4	55048	40	BR 5.0-6.0 (Blue W)	55013	3	Implassic FT3 3.25	55039	31
Osseospeed TX 3.0 S	55064	56	IR 3.5 (Yellow M)	55028	20	Implassic FT3 3.75	55049	41
Osseospeed TX 3.5-4.0 S	55046	38	IR 4.1-4.8 (Green R)	55028	20	Implassic FTP	55019	9
Osseospeed TX 4.5-5.0 S	55017	7	IR 4.8-5.5 (Blue W)	55028	20	Implassic HX3	55011	1
						Implassic TR2	55014	4
<b>Avinent</b>			<b>BlueSkyBio</b>			Implassic TRW2	55014	4
External, Coral & Ocean 3.5	55038	30	BIO Max NP Platform, All	55046	38			
External, Coral & Ocean 4.1	55011	1	BIO Max RP Platform 4.3-5.0	55037	29	<b>Dentalis</b>		
External, Coral 5.1	55013	3	BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	Regular	55051	43
Internal, Coral & Ocean 3.5, 4.1	55051	43	BIO Int Hex 3.5 Platform 4.7	55051	43			
Conical, Ocean 3.5	55066	58	BIO Int Hex 4.5 Platform 4.7-7.0	55051	43	<b>Dentegris</b>		
Conical, Ocean 4.1	55004	18	BIO One-stage Narrow, Regular, Wide	55014	4	3.3	55034	26
			BIO One-stage Short Regular & Wide	55014	4	3.75	55057	49
<b>B&amp;B Implant</b>			BIO Quattro Narrow 3.3	55061	53	4.5	55049	41
Duravit Slim	55021	11	BIO Quattro Regular 4.1 - 5.6	55062	54	3.8/4.1	55049	41
Duravit 3P	55035	27	BIO Conus 12, 3.0 Platform 3.0	55064	56	5.5	55020	10
Duravit Evolution	55035	27	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38			
			BIO Conus 12, 4.5 x 6mm	55046	38	<b>Dentis</b>		
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7	Cleanlant S-Clean 3.7 - 4.8	55020	10
			BIO Conus 12, 5.0 x 6mm	55046	38	S-Clean OneQ Narrow 3.3	55064	56
			BIO Trilobe 3.5	55016	6			
			BIO Trilobe 4.3-6.0	55017	7			

	No	Type		No	Type		No	Type
<b>Dentsply</b>			<b>Hahn</b>			<b>IRES</b>		
PrimeTaper 3.6	55046	38	Tapered 3.0	55058	50	External Hex Narrow 3.3	55012	2
PrimeTaper 4.2	55057	49	Tapered 3.5	55046	38	External Hex Regular 3.7 - 5.2	55011	1
PrimeTaper 4.8	55055	47	Tapered 4.3	55046	38	Internal Hex Narrow 3.3	55087	79
PrimeTaper 5.4	55055	47	Tapered 5.0	55037	29	Internal Hex Regular 3.7 - 5.2	55049	41
			Tapered 7.0	55048	40	Internal Octagon Regular 3.7 - 5.2	55014	4
<b>Dentium</b>			<b>Hi-Tec</b>			<b>J Dental Care</b>		
Superline	55037	29	Tapered Self Thread NTI 3.3, SLT 3,75	55034	26	JDEvolution 3.7-6.0	55067	59
NR Line	55066	58	Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	JDEvolution Plus+ 3.7-6.0	55067	59
Simple Line II GH 1,6	55014	4	Spark TPI	55034	26	JD Evolution S 3.25	55044	36
<b>DEVA</b>			Logic Plus	55034	26	JD ICON ultra, S 2.75	55064	56
Bone Level 3.5	55064	56	Logic LGI 3.0-3.5	55046	38	Conical Abutment	55080	72
Bone Level 4.5	55051	43	Logic LGI 4.3-5.0	55037	29			
<b>DIO</b>			<b>Hiossen</b>			<b>Keystone</b>		
SM Submerged 3.8	55049	41	HS system	55028	20	I-HEXMR 7,0-9,0	55051	43
SM Submerged 5.0	55020	10	ET system Mini D3.2	55047	39	Prima / Connex	55049	41
UF II Narrow	55064	56	ET system Mini D3.5	55015	5	Prima / Connex	55067	59
UF II Regular	55020	10	ET system Regular D4.0-7.0	55037	29	Prima / Connex	55041	33
			EK system Single D3.5	55046	38	Prima / Genesis Narrow	55049	41
			EK system Single D4.0-7.0	55044	36	Prima / Genesis Regular	55067	59
<b>Ditron Dental</b>			<b>IDI Evolution</b>			Prima / Genesis Wide	55041	33
Ultimate 3.75-6.0	55087	79	Standard implant 3.75-4.0	55004	18	Restore Narrow	55079	71
MPI 3.3-6.0	55087	79	B1One High	55080	72	Restore Regular	55011	1
			B1One Low	55033	25	Restore Wide	55013	3
<b>DSP Biomedical</b>			<b>IDI Implants</b>			TILOBEMAXX 7.0-9.0	55041	33
External Hexagon (HE) Standard	55014	4	BIO 3.5-5.0	55037	29			
Hexagonal Internal Standard (HIS)	55034	26	ID Cam ST 4.2	55010	13	<b>Klockner</b>		
Hexagonal Internal Morse (HIM)	55046	38	ID All	55010	13	VEGA Micro 3.0	55074	66
Hexagonal Internal Large (HiL)	55037	29	ID3 3.1	55021	11	VEGA Narrow 3.5	55065	57
Indexed Morse Taper (CMI)	55025	16	<b>Ids /Biotech dental</b>			VEGA Regular 4.0, 4.5	55034	26
			Kontakt, Kontakt-S 3.6 - 4.8	55045	37	KL Narrow 3.3	55012	2
<b>Easy Implant</b>			Reflect Aspire	55081	73	KL Regular 3.7, 4.2	55011	1
Ocean	55046	38	Reflect Certus	55019	9	KL Wide 4.7	55013	3
Lilas	55037	29	Reflect Tapered	55034	26	SK2	55038	30
			Reflect Recover	55020	10	Essential Cone	55067	59
<b>Edierre</b>			Reflect Rapid	55020	10	S3M 3.1	55079	71
3-lobe	55021	11	<b>Implance</b>					
4-lobe	55025	16	Tissue level	55014	4	<b>Konmet</b>		
			Bone level Regular	55037	29	Classic 4.0, 4.8	55014	4
<b>Euroteknika</b>			Bone level 3.3	55065	57	Cylindrical 3.5-4.0	55065	57
Internal Platform 3.0	55064	56	<b>Implant Direct</b>			Cylindrical 4.5-5.5	55057	49
Internal Platform 3.5	55081	73	Legacy 1/2/3/4 3.0/3.2	55012	2	Short Cylindrical 5.0-6.0	55078	70
Naturactis/Natural+ 3.0	55064	56	Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43	Conical 4.4/3.4	55065	57
Naturactis 3.5-5.0	55081	73	RePlant 3.5	55022	12	Conical 4.8/3.8, 5.5/4.5	55057	49
Uneva	55011	1	RePlant 4.3 - 5.0	55010	13			
			RePlant 6.0	55032	24	<b>Kyocera</b>		
<b>Futura implants</b>			ReActive 3.5 Platform	55022	12	FINESIA Bone Level 3.2-3.4 (NP)	55058	50
PRE 4.0	55014	4	ReActive 4.3 - 5.0 Platform	55010	13	FINESIA Bone Level 3.7-4.2 (RP)	55034	26
PRI 4.0	55037	29	ReActive 3.0 Platform	55046	38	FINESIA Bone Level 4.7-5.2 (WP)	55014	4
			Interactive 3.4 Platform	55037	29	FINESIA Tissue Level 3.7-4.7 (RP)	55014	4
<b>galimplant</b>			Screwindirect 5.0 PF 3.0 - 5.7	55070	62	FINESIA Tissue Level 4.7-5.2 (WP)	55014	4
iPx 4012	55046	38	<b>Implant One</b>			POI(HAC) EX 3.7 TP/ST	55085	77
			100 Series Std thread	55029	21	POI(HAC) EX 4.2 TP/ST	55025	16
<b>GC Implant</b>			200 Series Std thread	55029	21	POI(HAC) 2P 3.7/4.2	55004	18
Aadva Narrow	55081	73	300 Series Std thread (blue)	55029	21			
Aadva Regular	55066	58	400 Series Std thread (purple)	55035	27	<b>Lasak</b>		
Aadva Wide	55044	36	500 series Std thread (green)	55035	27	BioniQ QR	55066	58
Aadva short 4.2	55044	36	300 Series Wide thread (blue)	55029	21	BioniQ QN	55021	11
Genesis Plus 3.4-5.0	55085	77	400 Series Wide thread (purple)	55035	27			
SETIO Plus 3.8	55011	1	500 Series Wide thread (green) non-EU markets only	55129	95	<b>Leader</b>		
			<b>Implant Swiss</b>			Internal Hex, 3.3 - 5.5 mm	55034	26
<b>Geass</b>			Tissue level	55037	29	Cono-In	55067	59
Way Extra	55081	73	Bone level Regular	55065	57			
			Bone level 3.3			<b>Maxillent</b>		
<b>Global D</b>			<b>Implantvision</b>			iRaise	55049	41
In-Kone 4.0-5.0	55024	15	Internal Hex 3.7-5.7	55034	26			
TwinKon 3.5-4.5	55122	88	<b>Innobiosurg</b>			<b>Medentika</b>		
			MagiCore Regular	55018	8	Microcone NI	55021	11
<b>gmi</b>			<b>Intoss Anchor</b>			Microcone RI	55081	73
avantgard RP 3.75-4.25	55066	58	Internal Hex 3.7-5.7	55088	80	Quattrocone RI	55081	73
avantgard WP 4.75-5.75	55075	67				Quattrocone 30 AI	55081	73
frontier RP,WP	55051	43	<b>Intra-lock</b>			Procone 3.3	55119	85
phoenix NP 3.3	55038	30	Intrahex 3.4 - 4.5	55051	43	Procone 3.8	55031	23
phoenix Std 3.3 - Wide 5.0	55011	1	Gold & Blue 3.4 - 6.0	55015	5	Procone 4.3	55031	23
Monolith 3.0	55043	35	InDex 3.4 - 3.75	55020	10	Procone 5.0	55032	24
			InDex 4.0 - 4.3	55012	2			
<b>GP Implants</b>			FlatTop 3.75			<b>Medentis</b>		
SSI Spiral Shape implant	55051	43	<b>Innobiosurg</b>			ICX Premium/Active Master 3.75 - 4.8 (BL)	55062	54
			MagiCore Regular	55018	8	ICX Premium/Active Master 3.3 (BL)	55058	50
<b>gt medical</b>			<b>Intoss Anchor</b>			ICX Plus 3.45 Tissue Level	55039	31
HXE	55011	1	Internal Hex 3.7-5.7	55088	80			
HXE RP	55012	2				<b>Medical Instinct</b>		
HXE WP	55013	3				Hex 3.0	55071	63
HXI 4.1	55037	29				Hex 3.4, 4.0, 5.0	55052	44
002 HXI 4.7	55049	41				Cone 3.0, 3.4, 4.0	55071	63
OCT 4.3	55014	4				Cone 5.0	55072	64

	No	Type		No	Type		No	Type
<b>Medigma Biomedical</b>			<b>Nobel Biocare</b>			<b>PLATON JAPAN</b>		
Mars 3,5-6,0	55057	49	Brånemark System Ext hex RP	55011	1	Eight-Lobe Pro NC 3,5-4,1	55046	38
FixTite 3,5-6	55057	49	Brånemark System Ext hex NP	55012	2	Eight-Lobe Pro RC 4,6-5,0	55014	4
Fix-a-Dent 3,75-6	55057	49	Brånemark System Ext hex WP	55013	3	Type IV PRO	55020	10
			Tri-lobe Connection NP	55022	12			
<b>Medimecca</b>			Tri-lobe Connection RP, WP	55010	13	<b>P-I Brånemark</b>		
PT SLA Chaorum	55047	39	Tri-lobe Connection 6.0	55032	24	HEX 3.3	55115	81
			Conical Connection 3.0	55058	50	HEX 3.75 -5.0	55117	83
<b>Megagen</b>			Conical Connection Narrow	55046	38	AMP 3.3	55039	31
Anyridge	55051	43	Conical Connection Regular	55037	29	AMP 3.75 - 4.8	55118	84
Anyone Internal	55037	29	Conical Connection Wide 5.5	55048	40	MT 3.3 - 4.8	55046	38
Mini	55064	56	N1 Narrow	55127	93	Conical Abutment HEX, AMP, MT	55080	72
MUA	55080	72	N1 Regular	55126	92			
Octa	55121	87	MUA NP & RP	55080	72	<b>Prodent Italia</b>		
			On1 NP, RP, WP	55042	34	Twinner & Prime, Conomet, Collar: Orange	55027	19
<b>Meisinger Implants</b>						Twinner & Prime, Conomet, Collar:: Fuchsia	55051	43
OKTAGON Tissue Level RP, WP	55014	4	<b>Noris Medical</b>			Twinner & Prime, Conomet, Collar: Green	55067	59
OKTAGON Bone Level NC	55061	53	Tuff 3.3 - 6.0	55034	26	Twinner & Prime, Conomet, Collar: Yellow	55020	10
OKTAGON Bone Level RC	55062	54				Prime, Conomet, Collar: Blue	55020	10
			<b>Novem</b>					
<b>Microdent</b>			3.4	55039	31	<b>Qualia</b>		
Universal	55011	1	3.8	55016	6	NeoActive	55046	38
Microdent 2.8 - 3.25	55068	60	4.2, 5.0	55057	49			
Microdent 3.3 - 5.0	55011	1				<b>Quest</b>		
Microdent 5.5 - 8.0	55013	3	<b>NucleOss</b>			Hylock NP	55030	22
Genius 3.0	55025	16	Tpure Bone Level	55004	18	Hylock RP	55066	58
Genius 3.5 - 5.0	55051	43	T6 Bonel Level	55074	66			
Ektos 3.2	55023	14				<b>Resista</b>		
Ektos 3.7 - 8.0	55051	43	<b>OCO Biomedical</b>			WIS 3.0-5.0	55049	41
Trylogic 3.5 - 4.2	55051	43	Engage 3.25	55027	19			
Trylogic 5.0	55075	67	Engage 4.0-6.0	55051	43	<b>Rex Implants</b>		
			ERI 3.25	55069	61	External hex TL 1.8	55060	52
<b>MIS</b>			ERI 4.0-5.0	55014	4	Internal hex TL 2.9	55034	26
Conical connection Narrow (C1, V3)	55081	73	TSI 3.25	55068	60			
Conical connection Standard (C1, V3)	55087	79	TSI 4.0-5.0	55011	1	<b>Ritter</b>		
Conical connection Wide (C1)	55041	33				Multi Unit	55033	25
Seven and M4 Narrow	55027	19	<b>OsseoFuse</b>			Spiral SB/LA Narrow PF 2.9	55081	73
Seven and M4 Standard	55049	41	Hexaplus/ HexaPlus-S Sinus 3.75-7.5	55081	73	Spiral SB/LA Standard PF 3.75	55049	41
Seven Wide 5.0	55051	43	ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5			
Connect abutment NP, SP, WP	55115	81				<b>RS Implant</b>		
Multi Unit System	55033	25	<b>Osstem</b>			3.5 - 6.0	55020	10
			ETIII 3.5	55015	5			
<b>myplant GmbH</b>			ETIII 4.0-5.0	55037	29	<b>Schütz Dental</b>		
Myplant two 3,5-4,5	55077	69	US Regular	55011	1	Impla cylindrical	55052	44
Myplant Bio 4,0	55077	69	US Mini	55012	2	IMPLA micro Retention Cone 4.2	55052	44
			US Wide	55013	3			
<b>Neo Biotech</b>			US Wide PS	55011	1	<b>SGS Implant</b>		
IS II 4.0	55020	10	SS Regular, Wide	55028	20	P1, P7 Platform	55049	41
IS III 4.0	55020	10	TS Regular, Reg Short, Diam 6 & 7	55037	29	P1D, P5D, P7D	55049	41
Ektos 3.7 - 8.0	55051	43	TS Mini 3.5	55015	5	P7N	55059	51
Trylogic 3.5 - 4.2	55051	43	TS Mini 3.0	55047	39			
Trylogic 5.0	55075	67				<b>SHOFU</b>		
			<b>Osteo Plus</b>			BIOFIX 3.0	55064	56
<b>MIS</b>			Iridium PF 3.8 & 4.2	55079	71	BIOFIX 4.0	55046	38
Conical connection Narrow (C1, V3)	55081	73	She 3.45, 3.75	55011	1	BIOFIX 4.5	55046	38
Conical connection Standard (C1, V3)	55087	79	Shi 3.75, 4.2	55049	41			
Conical connection Wide (C1)	55041	33	Shi 5.0	55051	43	<b>SIC invent</b>		
Seven and M4 Narrow	55027	19				MUA	55120	86
Seven and M4 Standard	55049	41	<b>OT-Medical</b>			SICace 3.4 - 5.0	55052	44
Seven Wide 5.0	55051	43	OT-F2 3.4, 3.8	55053	45	SICmax, SICtapered 3.7 - 5.2	55052	44
Connect abutment NP, SP, WP	55115	81	OT-F2 4.1	55073	65	SICvantage® max, SICvantage tapered 3.0	55026	17
Multi Unit System	55033	25	OT-F2 5.0	55125	91	SICvantage® max, SICvantage tapered 3.7	55116	82
						SICvantage® max, SICvantage tapered 4.2-5.2	55054	46
<b>myplant GmbH</b>			<b>Overmed</b>					
Myplant two 3,5-4,5	55077	69	Easyline Easydip	55020	10	<b>Simply Implants</b>		
Myplant Bio 4,0	55077	69	ETLINE	55034	26	Internal Hex	55051	43
<b>Neo Biotech</b>			<b>Oxy Implant</b>			<b>SIN</b>		
IS II 4.0	55020	10	Internal	55020	10	Strong SW CM	55023	14
IS III 4.0	55020	10	Kone Conical Narrow	55021	11	Strong SW Int Hex	55035	27
			Kone Conical Regular	55051	43	Unitite CM 2.9	55084	76
<b>Neodent</b>			PSK Line Regular	55051	43	Unitite CM 3.5	55077	69
Cone Morse Connection	55025	16	External Mini	55079	71	Unitite CM Compact 4.0	55085	77
Grand Morse Connection	55066	58	External Regular	55011	1	Tryon 3.25 - 5.0	55011	1
Titamax WS Connection	55056	48	External Maxi	55011	1			
GM Micro Conical Abutment	55080	72				<b>Southern Implants</b>		
			<b>Paltop</b>			Tri-Nex 4.3-5.0	55010	13
<b>Neoss</b>			Standard 3.75-5.0	55034	26	Tri-Nex 3.5	55022	12
SP (∅ 3.5 implants and wider)	55004	18	Narrow	55027	19	External Hex	55011	1
NP (∅ 3.25)	55027	19	PCA Conical	55066	58	MAX Ext Hex	55011	1
Access Abutment	55060	52				Co-Axis Ext Hex 0°, 12°, 24°, 36°	55011	1
			<b>Phibo</b>			Piccolo Ext Hex	55050	42
			Aurea EVO NP	55027	19	IT-Connection	55014	4
			Aurea EVO RP, WP	55051	43	Deep Conical 3.0	55064	56
			TSA S3	55074	66	Deep Conical 3.5-4.0	55046	38
			TSA S4, S5	55034	26	Deep Conical Tapered 5.0	55017	7
			TSH S3, S4, S5	55011	1			

	No	Type		No	Type		No	Type
<b>Spiral Tech</b>			<b>Vulkan</b>					
ESi, Hex connection	55051	43	Internal Hex	55034	26			
Ultimate, Hex connection	55051	43	External Hex NP	55012	2			
ESi Conical 3.5	55044	36	External Hex RP	55011	1			
ESi Conical 4.3, 5.0	55037	29	Conical con MiP	55058	50			
			Conical con NP	55046	38			
<b>Straumann</b>			<b>Winsix</b>					
Tissue level, Regular & Wide Neck	55014	4	Internal Hex	55020	10			
Tissue level, Narrow Neck Crossfit (NNC)	55076	68						
Bone level / Bone Level Tapered Narrow (NC)	55061	53	<b>Xgate</b>					
Bone level / Bone Level Tapered Regular (RC)	55062	54	Conical 3.3 - 5.0	55034	26			
Bone level / Bone Level Tapered SC 2.9	55064	56						
Bone Level X (BLX) all platform	55027	19	<b>Xive</b>					
Screw retained Abutment NC, RC	55080	72	S-Plus 3.0	55030	22			
PURE Ceramic 4.1	55047	39	S-Plus 3.4	55061	53			
			S-Plus 3.8-5.5	55062	54			
<b>Surgikor</b>			TG Plus 3.4-4.5	55018	8			
Fixation 3.0	55058	50	<b>Zeramex</b>					
Fixation 3.5	55066	58	XT Implant 4.2, 5.5	55086	78			
Fixation 4.3	55037	29						
Solution 3.5, 4.0	55066	58	<b>Ziacom</b>					
Solution 4.5, 6.0	55037	29	Conical connection RP	55046	38			
Immediate, all	55049	41	Conical connection WP	55075	67			
Versatile, all	55049	41	External hexagonal connection NP	55079	71			
			External hexagonal connection RP	55117	83			
<b>Sweden &amp; Martina</b>			External hexagonal connection WP	55014	4			
Syra, Syra SP 3.9-5.0	55011	1	Internal hexagonal connection NP	55027	19			
Outlink <sup>2</sup> 3.75, 4.10, 5.0	55011	1	Internal hexagonal connection RP	55057	49			
Outlink2 3.3	55079	71	Internal hexagonal connection WP	55057	49			
Shelta 3.8-5.0	55025	16	<b>Zimmer</b>					
Premium Straight & SP 3.3-3.8	55025	16	Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41			
Kohno Straight & SP 3.3-3.8	55025	16	Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43			
Premium Straight & SP 4.25-6.0	55004	18	SwissPlus 3.8	55070	62			
Kohno Straight & SP 4.25-6.0	55004	18	SwissPlus 4.8	55014	4			
Prama 3.8 - 5.0	55079	71	Spline Twist 3.75-5.0	55038	30			
DAT Connection	55035	27	Eztetic NP 3.1	55065	57			
Shelta Abutment non-EU markets only	55128	94	<b>Zuga</b>					
			All	55027	19			
<b>Tag Dental</b>								
Axis	55051	43						
massif	55051	43						
<b>Tatum</b>								
Integrity 3.7 - 5.0	55015	5						
Integrity 6.0 - 8.0	55040	32						
T implant	55063	55						
<b>TBR</b>								
Infinity PF 3.5, PF 4.0	55004	18						
M Implant	55044	36						
<b>Thommen</b>								
SPI 3.5	55021	11						
SPI 4.0	55065	57						
SPI 4.5-6.0	55066	58						
<b>ticare</b>								
Osseus Mini, Standard	55011	1						
Osseus Maxi	55013	3						
Osseus Maxi ps	55011	1						
Inhex Mini	55021	11						
Inhex Standard	55046	38						
Inhex Maxi	55017	7						
Osseus Quattro Standard	55011	1						
Inhex Quattro Standard	55046	38						
<b>TRATE</b>								
ROOTFORM	55065	57						
<b>TRI Dental implants</b>								
Tissue level Octa	55014	4						
Tri-Vent, all	55034	26						
Narrow	55025	16						
<b>Trinon</b>								
Q2-Implant	55011	1						
QK-Implant	55014	4						
QZA 3.5 Platform	55049	41						
QZA 4.5, 5.7 Platform	55051	43						
<b>URIS</b>								
Omni Narrow 3.0-3.5	55115	81						
Omni Regular 4.0-6.5	55020	10						