



DNTLworks Equipment Corporation

Specifications

1310.19.00

ProSeal II Sealant Unit Lowspeed Handpiece 120V



| | |
|-------------------------|--------------------------------------|
| Amperage | 9.2 Amps |
| Compressor Type | Oil-Free Compressor/Vacuum Pump |
| Dimensions | 24.5"x10"x20 to ??" |
| Element Filtration | 50 micron Vacuum System |
| Flow Rate | 1 scfm @ 27psi Compressor |
| Flow Rate | 4.8 scfm @ 10 inhg Vacuum System |
| Hertz | 60hz |
| Horse Power | 1/16hp Compressor, 1/3hp Vacuum Pump |
| Operating Pressure | 65psi Air System |
| Shipping Size | 24"x24"x19" |
| Shipping Weight | 76 lbs |
| Sound Level | 46db @ 3' |
| Voltage | 120Vac |
| Volume Air Tank | 500ml |
| Volume Vacuum Container | 2qt, 1.9 Liter |
| Volume Water Container | 500ml |
| Wattage | 1100W to 1440W with Accessories |
| Weight | 60lbs |