



DNTLworks Equipment Corporation

Specifications

1500.07.00

ProAir I Oil Free Portable Air Compressor 120V



Amperage	8Amp
Compressor Type	Oil-Free
Flow Rate	2.75scfm @ 40psi
Hertz	60Hz
Horse Power	3/4Hp
Operating Pressure	70 to 100psi
Phase	1
Shipping Size	24"x24"x19"
Shipping Weight	58 lbs
Size	24-1/8"x 19-1/2"x8-3/4"
Sound Level	56db @ 3'
Voltage	120 Vac
Volume Air Tank	.95 US Gal.
Wattage	920W
Weight	42 Lbs