



DNTLworks Equipment Corporation Specifications

1540.11.00

ProStart Oil Free Portable Air Compressor 120V



Amperage	3.8Amp
Compressor Type	Oil-Free
Flow Rate	1.75scfm @ 40psi
Hertz	60Hz
Horse Power	1/3Hp
Operating Pressure	70 to 100psi
Phase	1
Shipping Size	26"x16"x19"
Shipping Weight	46lbs
Size	22"x14-1/2"x9-1/4"
Sound Level	52db @ 3'
Voltage	120 Vac
Volume Air Tank	.95 US Gal.
Wattage	437W
Weight	37 Lbs