



DNTLworks Equipment Corporation

Specifications

2200.34.00

ProCart II Self-Contained Treatment Console 120V



Amperage	10 Amps. 12 Amps Max with Accessories
Compressor Type	Compressor and Vacuum Pump: Oil Free
Dimensions	20"x21"x30"
Element Filtration	50 micron Vacuum System
Element Filtration	5 micron Air System
Flow Rate	2.5scfm @ 40psi. Compressor
Flow Rate	4.6scfm @ 10inhg. Vacuum
Font	Adobe Garamond Pro Regulator
Hertz	60hz
Horse Power	1/3Hp Compressor/Vacuum
Operating Pressure	70 to 100 psi Compressor
Regulated Pressure	65psi Air System
Regulated Pressure	40psi Water System
Shipping Size	23"x23"x38"
Shipping Weight	120Lbs
Sound Level	48db @3'
Voltage	120Vac
Volume Air Tank	2qt 1.9Liter
Volume Vacuum Container	2qt 1.9Liter
Volume Water Container	1 Liter
Wattage	1200 to 1440
Weight	96 LBS