

Prima Advance

Simplicity combined with the famous brand name provides effective passive mobilisation of the knee.

- Easy to use
- Tactile Adjustment Marks
- Available with Comfort Plastic Support
- Light
- Compact



Product Code: 4621008042

EASY TO USE

With simple and straightforward movement adjustments (angular limits, speed and patient sizing), the Kinetec Prima Advance is the most intuitive and quick to install device currently available.

It features tactile adjustment marks for patients with sight deficiency.

LIGHT

ONLY 12 Kg (24 pounds) to facilitate handling, the Kinetec Prima Advance is easy to transport making it ideal for the hospital and private physiotherapy environments.

COMPACTE

ONLY 94 cm (37 inches) in length to facilitate installation into a hospital bed or in the client's home and to aid storage.

RANGE OF MOTION From -5° to 115° of flexion. The Kinetec Prima Advance makes it possible to treat all pathologies of the knee.

Technical datas:

- Range of motion: -5° to 115°
- Speed (in °/minute): from 40°/min to 145°/min.
- Dimensions: length: 94cm, width: 33cm
- Sizes of the patient: from 1,45m to 1,95m
- Weight: 12kg
- Carriage width: 23cm
- Electrical power: 50VA
- Full leg: 71 to 99cm
- Tibia: 38 to 53cm, Femur: 33 to 46cm.
- Maximum patient body weight: 135 Kg.
- Simple hand control: Start/Stop/Reverse.
- The patient can stop and reverse the movement at any time.
- Voltage: from 100 to 240 Volts 50/60Hz.
- Type B classe I electrical device.

