


PMMA Resin Material

| Mono | Multi | Gum | Bi |



Mono Shade Temporary Restoration PMMA

CAD/CAM specialized temporary restoration material, composed of 99 polymethylmethacrylate (PMMA). With outstanding mechanical and compliance property, structural and shade stability, Mono PMMA can delivery authentic teeth like gloss , translucency, esthetically vivid restoration.



Indications:




Specification:

Bending Strength	≥120MPa	Dimension	12/14/16/18/20/22/25/30
System	Wieland/Amann/Zirkon/Sirona	Shade	Vita 16 shades/OM1/OM2/OM3
Water Sorption	≤25[μg/mm³]	Solubility	≤5[μg/mm³]

Multi-layer Temporary Restoration PMMA

CAD/CAM Anterior esthetic temporary restoration material, with revolutionary shade gradience technology, demarcation of shades realization, Multi-layer PMMA satisfies shade gradience from incisal to bottom. Excellent performance, stability, translucency and good bio-compatibility, multi-layer PMMA delivers no-stimulus features to the oral tissues.



Indication:




Specification:

Bending Strength	≥120MPa	Dimension	12/14/16/18/20/22/25/30
System	Wieland/Amann/Zirkon/Sirona	Shade	Vita 16 shades/OM1/OM2/OM3
Water Sorption	≤25[μg/mm³]	Solubility	≤5[μg/mm³]

Aidite / PMMA Resin Material

Gum Base Resin Material

With bio-fiber material, 100% imitate natural gum shades, effectively prolongs the denture durability. New technology also allows improvement on material bending, anti-impact and anti-fatigue performance.




Indications:

Bending Strength	≥120MPa	Dimension	12/14/16/18/20/22/25/30
System	Wieland/Amann/Zirkon/Sirona	Shade	Bionic Light Pink, Bionic Deep Pink, Bionic Yellow Pink, Bionic Brown Pink
Water sorption	≤25[μg/mm³]	Solubility	≤1.3[μg/mm³]

Full Contour Bi Shade Resin Material

Organically combined temporary PMMA with gum PMMA, Bi Shade PMMA is a resin material with excellent elasticity, strength, bio-compatibility and mechanical property.



Indications:

Bending Strength	≥120MPa	Dimension	20/22/25/30
System	Wieland/Amann/Zirkon/Sirona	Shade	Gum Base PMMA combination with Mono & Multi PMMA
Water sorption	≤25[μg/mm³]	Solubility	≤1.4[μg/mm³]



Bi-Shade PMMA



Bi-Shade PMMA