

One Drop



Two Tests



qLabs[®] Vet Coag Panel 2 PT/aPTT Combo

for Veterinarian Professional Use



- ❖ Bluetooth Communication
- ❖ aPTT Ratio
- ❖ Integral QR Code Reader

qLabs[®] Vet Coag Panel 2 PT/aPTT Combo

for Veterinarian Professional Use



REF number

QSV-1	Vet Coag Panel 2 test strips	12 strips/box
QSV-1-RC Pro	Re-Calcification Reagent (Accessory for Citrate whole blood testing)	12 vials/box

Feature Highlights

- Touchscreen handheld rapid testing platform with option of eStation II docking system
- Dry-reagent strips stored at room temperature for immediate use, no equilibration required
- Strong correlation to laboratory results
- Small blood sample volume, 10 μ L
- Stable external QC controls
- Built-in QR code scanner for accuracy and efficiency
- Built-in communication port for printers and computer networks
- Patented RACE[™] technology for accurate aPTT determination

Animal Species

- Canine (dog)
- Feline (cat)

Sample Type

- Fresh venous whole blood
- Citrate whole blood

Sample Volume

- 10 - 20 μ L

Time to Results

- 3 - 7 minutes

Reportable Time Range

- PT: 5 to 75 seconds
- aPTT: 6.5 to 500 seconds

Precision

- PT: In-house lot-to-lot and day-to-day precision \leq 5.0 CV% (mean value)
- aPTT: In-house lot-to-lot and day-to-day precision \leq 7.0 CV% (mean value)

Reference Intervals

- Canine (dog)
 - PT: 12 to 17 seconds
 - aPTT: 96 to 116 seconds
- Feline (cat)
 - PT: 13 to 20 seconds
 - aPTT: 96 to 124 seconds

Hematocrit Range

- 20% to 70%

Operating Temperature Range

- 10° C to 35° C

Operating Humidity Range

- 10% to 90% RH

Shelf Life at Room Temperature

- 12 months



CONTACT US

sales@ivetmedical.com
www.ivetmedical.com

Dubai

+971 50 523 8048
☎ +971 50 523 8048
DSO-IFZA, Dubai Silicon Oasis, UAE.

Lebanon

+961 1 497 740
☎ +961 81 310 308
3rd floor, Sukkariah Center, Dekwaneh, Lebanon.

Armenia

+374 93 720991
21, 1 Bagrevand St, Yerevan, Armenia.