



SpecializedRx Products, LLC  
7237 University Ave NE  
Fridley, MN. 55432 USA  
888.512.1209  
www.specializedrx.com

## SmartPopp<sup>™</sup> System Suggested Use Instructions

1. Separate mold cavities and visually inspect them, assuring they are free of any foreign particles.
2. SmartPopp were engineered to not require mold release lubricant when used with SpecializedRx Candy Bases. If needed, a very light coating of mold release (Pam<sup>®</sup> cooking spray, cocoa butter, etc.) can be applied to the interior of the cavities to help aid lollipop removal.
3. Insert sticks provided into the SmartPopp<sup>™</sup> molds. Be sure to use only Specialized Rx provided sticks as other sticks should not be used due to dimensional/quality requirements to function with molds.
4. Place SmartPopp<sup>™</sup> Stand on a clean, level working surface.
5. Then insert molds into the preferred open spaces of the stand. It should be a tight fit as this is needed to maintain pressure to both halves of the SmartPopp<sup>™</sup> while assembled for filling.
6. Then visually inspect to assure each stick is resting on the level working surface. This will ensure consistent filling accuracy from unit to unit. If a stick is not resting on the working surface, pull the stick down to meet the surface.
7. Once all molds and sticks are in place, they are ready to start filling. We recommend using the variety of Specialized Rx bases available found online at [www.specializedrx.com](http://www.specializedrx.com)
8. Carefully fill each mold to prevent overflow. For ease of filling, recommended use of a catheter tip style syringe (Item E0568) or beaker based on preference.
9. Once mold contents are cooled and set, the molds can be removed from the stand by pushing each unit up from the bottom. Lollipops can be inspected by separating the halves to ensure contents move freely. Individual SmartPopp<sup>™</sup> can be separated by gently breaking the two center tabs off each side.
10. Each unit should be checked for proper weight according to your quality control procedure to maintain appropriate weight tolerance.
11. Once the weight has been verified, the child resistant caps can be applied to each unit. Use only Specialized Rx supplied child resistant caps. Other caps should not be used due to dimensional requirements to securely fasten caps to SmartPopp<sup>™</sup> molds.

### Note:

- Procedures should be adjusted base on use and the variety / quantity of ingredients to be used.
- Use only Specialized Rx supplied lollipop sticks and child resistant caps. Use of any other product has not been verified for functional use or child resistant assurance. Only Specialized Rx supplied SmartPopp<sup>™</sup> System components meet the requirements of the CFR and ISO standards for child resistance.

**Disclaimer:** Specialized Rx makes no warranties or claims to the accuracy or completeness of the above information given the variations in use and formulations. These suggested use instructions are provided as user guidance considering varying circumstances. The SmartPopp<sup>™</sup> System was designed for compounding single use units. Components should be discarded or recycled following patient use.

### Associated Products:

SmartPopp CR Lollipop Molds, 24 Count, Item# E0500-01  
SmartPopp Prescription Lollipop Stand (30 count) Item# E0502-00  
Syringe, Sterile Catheter Tip (60ml) Item# E0568-00 (1pk) / E0568-03 (5pk)  
Nextol SF<sup>™</sup> Candy Base w/ Bitter-Bloc (sugar-free) 1Kg, Item# B0030-55  
Nextol SF<sup>™</sup> Candy Base w/ Bitter-Bloc (sugar-free), 2.0Kg, Item# B0030-62  
Krystal SF<sup>™</sup> Quickset Candy Base w/ Bitter-Bloc (sugar-free)1Kg, Item# B0020-55  
Krystal SF<sup>™</sup> Quickset Candy Base w/ Bitter-Bloc (sugar-free)2.0Kg, Item# B0020-62