

Tablet Triturate Molds, 100 Cavity Hard Coated

Two-piece 100 cavity molds are designed for making tablet triturates or troches. Each mold consists of a base plate with alignment pegs and a top plate with corresponding holes. To make tablets, place the top plate on a Glass Pill Tile or Ointment Slab then fill the holes evenly with the tablet paste using a hard rubber spatula to compress the mixture. When the paste has hardened, slide mold off edge of slab, align the top plate alignment holes (left and right side) over the base plate alignment pegs (left and right side) and press down to punch out tablets from top plate. Allow to dry/harden.

Made of durable, corrosion and wear-resistant **hard-coated** aluminum. Hard-coat is far superior to other "anodized" molds on the market in terms of scratch resistance, durability, and general use surface wear. Available in three sizes: 60, 100 and 200 mg. All are 6 x 4.5 in. / 15.2 x 11.4 cm. **Patent pending.**

Mold Size	60 mg	100 mg	200 mg
Hole Size \varnothing x D:	0.18 x 0.13 in. / 4.6 x 3.2 mm	0.24 x 0.13 in. / 6.0 x 3.2 mm	0.24 x 0.25 in. / 6.0 x 6.4 mm
Hole Volume:	0.053 ml	0.090 ml	0.181 ml

Cleaning: Separate all plates before cleaning. Recommend hand cleaning using a mild detergent – Recommend using **Relianox™ Critical Cleaning Liquid Detergent (SKU# P3023)**.

Relianox™ is the ideal critical cleaning liquid detergent for use in laboratories, clinics, hospitals, pharmacies, veterinary, cleanrooms, and industrial facilities. Relianox™ is a safe commercial grade detergent that is effective in cleaning glassware, instruments, rubber, plastic, and metal components, surfaces, equipment, tanks, and numerous other applications. Relianox™ is a phosphate and soap-free concentrated detergent that is easy to use and offers superior activity for highly effective cleaning in critical environments.

Cleaning Agents: Recommend **Relianox™ Critical Cleaning Liquid Detergent (SKU# P3023)**. Do not use acidic or caustic cleaning solutions. Do not clean with abrasive cloth, gritty materials, or steel wool. Use only plastic brushes or nylon/non-abrasive pads. Hard-coated aluminum may be cleaned in an automatic dishwasher using only a mild detergent.

Isopropyl Alcohol Cleaning or Drying: May be cleaned with isopropyl alcohol. Wipe with a cloth dipped in isopropyl alcohol, then immediately wipe dry with a clean cloth. May be dried in a hot air oven or using a portable hair dryer.

Autoclave safe.



SpecializedRx
 3787 95th Ave NE, STE 100
 Circle Pines, MN. 55014 USA
 888.512.1209
 www.specializedrx.com

TritTab™ SF/LF (Sugar-Free/Lactose-Free) Tablet Triturate Base

CATEGORY	TritTab SF/LF was specifically engineered for molded tablet triturates with the ability to use alcohol alone as the wetting solution and quickly disintegrate in the presence of moisture. TritTab SF/LF is formulated to be an extremely versatile tablet triturate base that is sugar-free, lactose-free, exhibits minimal contraction, excellent friability properties, and compatible with many drugs that react chemically with sugars. This base contains only natural ingredients, sweetened with our TruClear Stevia Plus, and contains our natural Bitter-Bloc™ Technology for optimal palatability. Bulk Density, at 60°C approx. 0.47 g/cm ³ (measured).				
PRICING	SKU	Size	Price	NDC	**Inquire for larger sizes
	B0215-40	100g	\$20.00	72467-0215-40	
	B0215-50	500g	\$45.00	72467-0215-50	
	B0215-62	2Kg	\$125.00	72467-0215-62	
	B0215-65	5Kg**	\$265.00	72467-0215-65	

TritTab™ Velvet / Sugar-Free/ Fast Dissolving Tablet Triturate Base

CATEGORY	TritTab Velvet was formulated as a convenient and versatile sugar-free tablet triturate base utilizing specifically sourced ingredients and particle sizes for compounding a wide range of APIs appropriate for this dosage form. Specifically engineered for molded tablet triturates requiring an alcohol/water solution as the wetting agent that will quickly disintegrate when placed in the mouth or buccal cavity. TritTab Velvet exhibits optimal flow characteristics, minimal contraction upon drying and excellent friability properties in tablet form. This base includes our natural Bitter-Bloc™ Technology for optimal palatability. Bulk Density, at 60°C approx. 0.64 g/cm ³ (measured).				
PRICING	SKU	Size	Price	NDC	**Inquire for larger sizes
	B0216-40	100g	\$17.00	72467-0216-40	
	B0216-50	500g	\$42.00	72467-0216-50	
	B0216-62	2Kg	\$115.00	72467-0216-62	
	B0216-65	5Kg**	\$255.00	72467-0216-65	

TritTab™ Natural Liquid Flavor Concentrate

CATEGORY	Formulated with natural extracts for use in either TritTab™ SF/LF or TritTab™ Velvet as the wetting agent. Can be used alone as the wetting agent with TritTab™ SF/LF tablet triturate base. Replaces the high-proof alcohol when compounding tablet triturates with either TritTab™ Bases. Water-soluble. 24-month BUD from date of manufacture. The hazmat shipping fee applies to this product. See current shipping policy online. Some shipping restrictions apply.			
PRICING	TritTab™ Peppermint Ice			
	SKU	Size	Price	NDC
	F0302-80	500 ml	\$38.00	72467-0302-80
	F0302-85	3.8 L	\$130.00	72467-0302-85