# takes you further



# Ergonomically and stylishly designed rollators

Living life to the fullest, that is what Rollz is all about. Be it for a walk in the park, a visit to a museum, a holiday trip or a shopping afternoon. Our ergonomically and stylishly designed Rollz Motion<sup>2</sup> and Rollz Flex will bring you further, while also supporting an active and independent life.





## Enjoy adventures with Rollz Motion<sup>2</sup>

## The 2-in-1 rollator and wheelchair

The Rollz Motion<sup>2</sup> is perfect for those who want to continue their journey but worry about how far they can go. This 2-in-1 mobility equipment allows you to walk in comfort and also benefit from the ability to transform it into a wheelchair when needed.



## Product features Rollz Motion<sup>2</sup>

## **Stylish**

Award winning Dutch design rollator and wheelchair in one, for going out in style.

## Safe

Innovative drum brakes, suitable to wet and dry conditions.

## Convenient

Rest on a comfortable cushioned seat or transform the rollator into a wheelchair.

## Easy aid

Cross doorsteps and kerbs easily with the kerb assistant.



## Maneuverable

Light steering and easy motion with the maneuverable front wheels.

## Stable

Walk tall and steadily between the ergonomic hand grips that are adjustable in height.

## **Compact**

The Rollz Motion<sup>2</sup> is easily foldable into a compact package that fits in almost any car.

## Comfortable

The large rear wheels with soft tyres and foam handles ensure a comfortable walk.








## Easy to transform

With a flick of the wrist you can transform the Rollz Motion<sup>2</sup> rollator into a wheelchair and back again. This will give you confidence to go further and rest along the way, while enjoying a day out with your family or friends.

The **Rollz Motion**<sup>2</sup> is available in four colours: pebble white, dark purple, matt black and island blue.

The **Rollz Motion**<sup>2</sup> **Small** is available in island blue or matt black.









## Smooth and safe ride on uneven terrain

## All terrain rollator and wheelchair in one

Thanks to the four pneumatic tyres, the Rollz Motion Performance provides a comfortable ride, even on bumpy roads. The extra thick padding of the armrests and the ergonomically shaped handles add to the comfort. Moreover, the footrests are adjustable in height.



## Product features Rollz Motion Performance

## Stylish 2 in 1

Dutch design for going out in style, which transforms from rollator to wheelchair within seconds.

## **Ergonomic grips**

Ergonomically shaped soft handles for optimal grip.

## **Eye-catching color**

Modern matt coating in the new color jungle green.

## Visible

Visible in the dark thanks to multiple reflective parts.



## Soft padding

Armrests with extra thick foam padding to increase comfort.

## Firm seat

Wedge-shaped cushion for a better sitting position.

## Air filled tyres

Sturdy pneumatic tyres enable more comfortable rides on bumpy roads.

## **Stable**

Walk tall and steadily, while safely supported by a sturdy design.






## **Specifications**

	Rollz Motion <sup>2</sup>	Rollz Motion <sup>2</sup> Small	Rollz Motion Performance
Rollator weight	11,6 kg	11,6 kg	11,4 kg
Wheelchair weight	3,4 kg	3,4 kg	3,4 kg
Maximum weight capacity	125 kg	125 kg	125 kg
User height	1,60 - 1,90 m	1,50 - 1,85 m	1,60 - 1,90 m
Seat height	55 cm	50 cm	55 cm
Handle height range	87-98 cm	82-94 cm	86-98 cm
Seat size rollator	20 x 45 cm	20 x 45 cm	20 x 45 cm
Seat size wheelchair	42 x 45 cm	42 x 45 cm	42 x 45 cm
Overall width	67 cm	67 cm	67 cm
Overall length	63 cm	63 cm	63 cm
Width folded	31 cm	31 cm	29 cm
Front wheel diameter	8 inch	8 inch	Tyre 8 x 1¼ inch
Rear wheel diameter	12 inch	12 inch	Tyre 12½ x 2¼ inch

Warranty: 5 years on the frame and 2 years on other parts with exception of normal wear and tear.





# On the go with Rollz Flex

## Design rollator with large bag

At first glance, it does not appear to be a walking aid. This was exactly the intention of the designers, to create an enjoyable design for going outdoors. This rollator maneuvers easily and you can walk stable and comfortable, while the bag carries your groceries.





## Maneuverable

You can easily maneuver anywhere, even in between obstacles. The front and rear wheels are set close to each other, and this way your rollator rotates close around you. It allows you to turn smoothly in a narrow space.

If you make a turn, the wheels do not get in the way. You can step aside easily and take something out of a closet or rack. The Rollz Flex functions the same even when the bag is loaded to its capacity.

With ease you can go up the kerbs with the help of the kerb assistant.





## Stylish

## Product features Rollz Flex

Fully certified rollator with the look of a designer shopper.

## Flexible

Easy adjustable handlebar for different heights.

## **Practical**

The large water-resistant bag can carry up to 20 kg. Includes removable inner bag.

## Comfortable

Comfortable ride thanks to the soft tyres.

## Maneuverable

Light steering and easy motion thanks to the small turning cycle.



## Safe

Safe drum braking in any weather conditions. The brake cables are concealed within the frame.

## Convenient

Rest on the leather look seat while having the handlebar as a supportive backrest.


## Easy aid

Go over doorsteps and kerbs easily with the kerb assistant, even with a full loaded bag.

## **Light and compact**

The Rollz Flex weighs only 7.4 kg and is easily foldable into a compact package.





## **Product details**

	Rollz Flex	Rollz Flex Small
Rollator weight	7,4 kg	6,8 kg
User height	1,55-1,95m	1,40-1,60 m
Seat height	62 cm	52.5 cm
Seat width	44 cm	44 cm
Weight capacity	125 kg	125 kg
Overall width	62,3 cm	62,3 cm
Overall length	60 cm	60 cm
Depth folded	33 cm	33 cm
Height	80-100 cm	69-81 cm
Bag capacity	25 litre/20 kg	14 litre/20 kg
Brake type	classic, flip	classic
Wheel diameter	8 inch	8 inch

The **Rollz Flex** is available in two frame colours: matt black and white. The bag is made out of sustainable and water resistant material and is available in two colours (bright purple and denim grey).

Warranty: 5 years on the frame and 2 years on other parts with exception of normal wear and tear.

In rollator position the **Rollz Flex** is a fully TÜV certified rollator.











Due to the unique design of the Rollz travel cover, both the Rollz Motion and the wheelchair package are incorporated into one compact package. Using the travel cover will help preventing damages.

## **Rollz Motion shopper**

This stylish and functional shopping bag can be clicked onto the Rollz Motion frame. You can put your shoppings directly in the shopper and take it along, so you won't have to carry a shopping bag.





The Rollz cup holder is very easy to use. Thanks to the sliding system, the holder can be placed into the desired position quickly, enabling the user to take a water bottle along.



## Rollz Motion 3-in-1 wheelchair package holder

With these holders you can also hang your wheelchair package over the seat. This gives you more walking space and freedom of movement. A cane and shopping bag holder are integrated.



## Rollz Motion umbrella

This umbrella protects you from rain and sunshine. It gets attached to the frame so you can use both hands on the rollator for best support.

## **Rollz Motion back support**

Out on your own and left the wheelchair package at home? This back support offers you a short rest, while leaning against the soft strap.



### **Rollz Motion seat belt**

The seat belt keeps you fixed in the wheelchair for extra safety. It can be easily placed on the wheelchair package and adjusted to the desired length.



## Rollz Flex cane holder

The cane holder attaches to the handlebar and stays there. The corresponding part has been designed to fit almost any walking stick.



## **Rollz Flex shopping hooks**

If you go shopping, this set of hooks offers the possibility to carry along another shopping bag.



### Rollz Flex umbrella

This umbrella attaches to the frame and protects you from rain and sun. You can still use both hands on the walker for support.



## Rollz International BV www.rollz.com













zint

RollzInternationa