

TECH WEST ECOVAC DRY VACUUM PUMP

Our EcoVac system provides the advantages of proven technology, durable construction, quiet operation and oil-free performance with the latest in water and energy-saving features. This is not only important to you as an environmentally responsible business owner, it also prepares your practice for the upcoming change to metered water usage in many areas of the country. We utilize a Variable Frequency Drive (VFD) to control the vacuum levels. With no suction being used, the motor is running one third the speed and using very little power. Every time a suction hose is being used the motor automatically speeds up to maintain constant suction.

Maintenance free unit with no traps to clean, no oil to check or change. Truly a ZERO maintenance machine.



ECO DRY-VAC SYSTEMS

MODEL NUMBER	MAX USERS	# OF HEADS	TOTAL H.P.	CFM @ 8 HG	TANK GALLONS	VOLTAGE	TOTAL AMPS	BREAKER SIZE	DB(A)	EQUIP. WEIGHT	EQUIP. HEIGHT	EQUIP. WIDTH	EQUIP. DEPTH	SHIPPING WEIGHT	SHIPPING HEIGHT	SHIPPING WIDTH	SHIPPING DEPTH	LIST PRICE
VPD2S2	2-3	1	1	20	14	208/230	10	20	60	139	52	23	20					
VPD4D2	4-6	2	2	40	14	208/230	10	2X20	62	282	39	40	20					
VPD5S2	5-7	1	2	45	14	208/230	12	20	62	139	52	23	20					

DVU-S Upgrade 2-3 user single to 5-7 user: 1,855.00

DVU-D Upgrade 4-6 user dual to 10-14 user: 3,710.00

FDA
RECOGNIZED

5 YEAR WARRANTY OR 5000 HRS WHICH EVER COMES FIRST

Meets or exceeds NFPA applicable guidelines