



cara i700 wireless*

The power of intraoral scanning is now even more convenient. Convince yourself and start working wirelessly today!

* powered by Medit.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

cara i700 wireless*

WIRELESS No wire, no hassle	STABLE & SPEEDY Up to 70 FPS	FOR LONGER USE Battery life of up to 8 hours	
COMPACT & LIGHT 328 g	REMOTE CONTROL Mode	UV-C LED Disinfection	180° Reversible tip

Category		cara i700 wireless
Scanning Technology	Scanning Frame	Up to 70 FPS
	Imaging Technology	3D-in-motion video technology/3D full color streaming capture
	LED Source	LED
	Anti-fogging Technology	Adaptive anti-fogging
Accuracy	Full-arch	10.9µm ± 0.98
Handpiece	Dimensions	313 x 44 x 47.4 mm
	Total Weight	328g (Includes battery, tip)
Tip	Tip Size	22.5x17.1 mm
	Scan Area	15x13 mm
	Autoclavable	Autoclave up to 150 times Autoclave 121°C 30 mins Autoclave 134°C 4 mins (Pre-Vacuum Type) Autoclave 135°C 10 mins (Gravity Type)
	Reversible Tip	Yes
Special Mode	Remote Control Mode	Yes
	UV-C LED Disinfection	Yes
Connectivity		Wireless USB 3.1 Gen 1 (Wireless Hub)

- Performance can vary depending on PC specs
- Accuracy test of i700 wireless scanner, Medit R&D test

* powered by Medit.