


NEW TECH WEST

ELITE SERIES ULTRA CLEAN OILLESS COMPRESSOR

Our Elite Series Oilless Air Compressor system provides the advantages of proven technology, durable construction, quiet operation and oil-free performance with the latest in energy-saving features. This is important to you as an environmentally responsible business owner.



ELITE SERIES OILLESS AIR COMPRESSOR SERIES

MODEL NUMBER	MAX USERS	AIR DRYER	QUIET COVER	# OF HEADS	TOTAL H.P.	CFM	TANK GALLONS	VOLTAGE	TOTAL AMPS	BREAKER SIZE	DB(A)	EQUIP. WEIGHT	EQUIP. HEIGHT	EQUIP. WIDTH	EQUIP. DEPTH	SHIPPING WEIGHT	LIST PRICES
ACO2S2	2-3	YES	NO	1	1.5	5.5	10	208/230	8	20	70	180	30	28	20.5		
ESCS-100	QUIET COVER: OPTIONAL FOR THE ACO2S2										68	50	32	28	20.5		
ACO4D2	5-7	YES	NO	2	3	11	20	208/230	16	20	72	248	44	42	21		
ACO6T2	8-10	YES	NO	3	4.5	16.5	30	208/230	25.5	40	74	320	46	41	31		
ESCD-100	QUIET COVER: OPTIONAL FOR THE ACO4D2										68	50	38	29	33		

A 5 YEAR or 3500 HR MOTOR WARRANTY. All other components carry a 2 year warranty
Meets or exceeds NFPA applicable guidelines

FDA
RECOGNIZED