
REACHER GRABBER TOOL

User Manual



Review the user manual carefully prior to first use. To ensure proper product operation, pay special attention to the safety precautions. Retain this user manual for further reference.

SAFETY PRECAUTIONS

To avoid damage or injury to the user or third parties, read all instructions and safety recommendations prior to use of this product. Continue to follow all instructions and guidelines.

WARNING:

- Keep out of reach of children.
- Do not leave children or minors unattended while using this product.
- Do not use it for other means than those described in this manual.

PRODUCT SPECIFICATIONS:

Length:	32.3" (82cm)
Weight:	0.25 lb (114g)
Opening:	4" (10 cm)
Max. Load Capacity:	4 lbs (2 kg)
Max. Temperature:	140° F (60 °C)

PRODUCT INFORMATION & SUGGESTED USE:

- Ideal for grabbing objects in hard-to-reach places such as: floor, ground, behind furniture, under beds or high shelves.
- The 4-inch jaw opening has interlocking tips for a secure grip. It's easy to pull when closed, even for people with low hand strength.
- The tip contains a convenient magnet which you can use to pick up light metal objects.
- Can be easily pulled apart in the middle and clipped together for easy portability.
- Includes plastic clip to attach to walkers, wheelchairs and canes.

CLEANING AND MAINTENANCE:

All parts can be wiped clean using a cloth and a gentle household detergent or disinfectant.

LIMITED LIFETIME WARRANTY:

Your Medic Therapeutics product is backed by a limited lifetime manufacturer's warranty. Medic Therapeutics will repair or replace your device at any time should it fail due to a defect in material or workmanship, subject to the certain limitations.

This limited warranty does not cover any damage that results from unauthorized or improper use, service, or repair. Further, it does not cover damage caused by accident, impact, negligence, or normal wear and tear. Should you discover your Medic Therapeutics product is not functioning properly, please send your device to our repair center for evaluation. If your product cannot be repaired or serviced, we will reserve the right to change it for a similar or newer model.

Please note that a flat rate of \$35.00 will be charged to cover a service evaluation fee and return shipping of your device. All warranty claims must be accompanied by a copy of your proof of purchase from an authorized retailer. Please send your device, proof of purchase, and a check or money order in the amount of \$35.00 made out to Medic Therapeutics to:

Address:

Medic Therapeutics Service Center
3069 Taft Street
Hollywood, FL 33009

Contact:

warranty@medictherapeutics.com