



Email: admin@synemed.info



# Capno Fit **ETCO2 Sampling Products**





#### NINGBO SHENGYURUI MEDICAL APPLIANCES CO.,LTD.

Company Address: No. 138, Binhaisi Rd., Hangzhou Bay, Cixi, Ningbo, Zhejiang, China Tel: 86-574-63006301-8015 Fax: 86-574-63259048 Mobile: 86-13586615067 E-mail: sales@soundway-medical.com www.soundway-medical.com

## Capno Fit ETCO<sub>2</sub> Mask

ETCO, mask is a device ideal for patient convert to mouth breathing.

- Incorporate a bevel tip in the dome of mask Ensure exhale CO, gas sampling with simultaneous insufflation of O<sub>2</sub> or gaseous analgesia
- Co-extruded connecting tubing Avoid interfere during patient treatment if the tubing is twisted

- Transparent and clear mask Ensure clear observation of patient for accurate visualization of cyanosis
- · Optional hydrophobic filter



#### **Order Information**





Dank Normalian	Description	Tubing Length	Luer Connector		Corr. Ot.
Part Number			Male	Female	Case Qty
C0301	ETCO <sub>2</sub> Mask, Large Adult	7FT	•		50
C0302	ETCO <sub>2</sub> Mask, Small Adult	7FT	•		50
C0303	ETCO <sub>2</sub> Mask, Pediatric	7FT			50



## Sampling Line

- Co-Extruded PVC/PE sampling lines Ensure more rapid and accurate
- · Optional hydrophobic filter

- · Choice of male or female luer connector for machine fitting
- Standard PVC sampling lines available

### **Order Information**





Part Number	Description	<b>+</b> 1	Luer Connector		6 0
		Tubing Length	Male	Female	Case Qty
C0324	CO <sub>2</sub> Sampling Line	7FT	•	•	100
C0325	CO <sub>2</sub> Sampling Line	10FT	•	•	100
C0326F	CO, Sampling Line	10FT	2•		100