

Retract Micro Dissection Needle

Retract is a small scale fingerswitch that contains an ultra sharp tungsten needle. The considered ergonomics of the device offering a greater level of accuracy in microsurgical procedures, especially when working around the eyes.

The soft touch buttons are placed closer to the operative site to offer optimal control, more like writing with a pen.

Retract is furthermore a safety device - the needle lives inside the handle and is only extended to use the device, then can be retracted when not in use and for safe disposal.

The Surgeon's 'ballpen'

Control is absolute, so close to the operating tip that it becomes an extension of the finger, no bending or lateral movements, perfect balance between the user and microdissection surgery.


The Nurse's friend

Needlestick injuries are minimised with the retracting needle: press to release and twist to retract back into the safety of the handle

Three products in one

Two needles and handle all in one product

RETRACT
micro dissection needle



Product REF	QTY
PMS8001 Retract electrosurgery fingerswitch with retractable micro dissection needle	10