Reynard IPA Surface Disinfection Wipes



IPA Surface Disinfection Wipes

Reynard's IPA Surface Disinfection Wipes are ideal for all healthcare professionals who require an effective solution to their infection control needs. The large sized wipes enable a quick and efficient disinfection to a variety of washable hard surfaces such as tables, counter tops, trolleys and stainless steel worktops.



Directions: Use the required number of wipes. Dispose of used wipes properly; do not flush.



Disinfection

Hygienic, easy to use, convenient & efficient. Reynard IPA Surface Disinfection wipes will disinfect most hard surfaces and one easy step.

Fast and effective

Reynard Surface Disinfection Wipes are impregnated with 70% Isopropanol and a very low level of non-ionic surfactant which provides effective disinfection of clinical surfaces

Cost-effective

Reynard IPA Surface Disinfection wipes will disinfect in one step. They are highly portable, convenient and practical.

Technical Performance:

This product has been tested and proved to be effective against the following micro-organisms:

- MRSA
- E.Coli
- Salmonella Typhimurium
- Listeria moncytogenes
- P.aeruginosa
- Proteus vulgaris
- C. albicans
- A. niger

These IPA Surface Disinfection products have been tested as per EU Standards under "Dirty Conditions" to EN1276

Reynard's IPA Wipes have NZSFA (MPI) C44 Approval

Product specifications and availability may vary according to country.

Always read the label and product information before use.

Supporting documentation: Full technical specifications are available upon request. **Bespoke solutions:** Should your market require more or fewer wipes, lightly scented alternatives or even biodegradable wipes, for example. Let us know. We are more than happy to develop bespoke solutions.

IPA Disinfectant Wipes



One step disinfecting. Reynard IPA Surface Disinfection Wipes are impregnated with 70% Isopropanol Presented in the familiar canister format they are highly portable, convenient and practical.

The flip top lids prevents these wipes from drying out.



RHS 402 IPA Surface Disinfection Wipes

Technical specifications

- 22 x 20cm, 33 gsm non-woven
- 160 wipes per canister,
 12 canisters per carton

Dispenser

• Perspex wall mounted bracket



The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. Regarding material compatibility are only possible in particular cases. Our recommendations are without obligation and do not constitute a warranty. They do not preclude the company's own testing for the intended purposes and processes. In this respect we cannot accept any liability.

This complies with our general conditions of sale and supply.

United Kingdom

Boston House
Grove Technology Park
Downsview Road, Wantage
Oxfordshire OX12 9FF
+44 845 263 7335

Australia

Suite 403, 44 Hampden Road Artarmon Sydney, NSW 2064 enquiries @reynardhealth.com 1300 667 199

New Zealand

17 Napier Road Havelock North Hawke's Bay 4130 New Zealand +64 (0)6 650 0709



Tel:(852) 3105 2821 Fax:(852) 3105 0197 E-mail: admin@synemed.info