


# ProSling

## Patient lifting slings for **yoke** frame lifters



A range of slings designed to lift patients from seated (chair) or reclined (bed or floor) positions.



- Versatile enough for general patient transfers, toileting or hygiene care
- Manufactured from the highest quality material; conforms to industry standards (AS/NZS ISO 10535:2011)
- Heavy duty, quick dry 3 mm quilted foam padding in leg and head support sections delivering increased patient comfort
- Brightly coloured handles and adjusters, for easy access and visibility
- Colour coded sling straps provide easy identification of sling size
- Three loop support strap system, controls positioning of head and leg supports
- Suits users to a maximum weight of 205 kg (Small and Medium); and 300 kg (Large, Extra Large and Extra Extra Large)

Will suit the majority of lifter products on the market. Some lift manufacturers make claims that only their brand of sling should be used on their lifter. There is no Australian Government or TGA directive to support these claims.



## General Purpose



Provides good support for patients with basic upper body control.

## General Purpose with Head Support



In built head support creates a secure transfer for patients with less controlled body movements.

## Hygiene




Designed for toilet transfer, includes a waist strap with buckle fastening for added security and generous padding under arms and around mid section.

## Hygiene with Head Support



Designed for toilet transfer, in built head support creates a secure transfer for patients with less controlled body movements.



### DIMENSIONS in mm

## General Purpose



SIZE	A WIDTH	B LEG WIDTH	C LEG LENGTH	D DEPTH	E LENGTH	F HEAD STRAP	G LEG STRAP	H TORSO BAND	CAPACITY	CODE
Small	840	230	670	490	1160	450	450	n/a	205 kg	SL147181
Medium	890	240	655	545	1200	450	450	n/a	205 kg	SL147281
Large	1135	275	685	640	1325	450	440	n/a	300 kg	SL147381

## General Purpose with Head Support



Small	830	230	540	790	1070	300	450	n/a	205 kg	SL146181
Medium	950	240	550	920	1220	300	450	n/a	205 kg	SL146281
Large	1155	260	700	1000	1350	300	450	n/a	300 kg	SL146381
Extra Large	1350	260	750	1150	1500	300	450	n/a	300 kg	SL146481
Extra Extra Large	1450	270	800	1200	1600	300	450	n/a	300 kg	SL146581

## Hygiene



Small	940	175	625	195	820	320	645	890	205 kg	SL267181
Medium	1040	175	695	195	890	320	645	990	205 kg	SL267281
Large	1170	175	740	195	935	320	655	1120	300 kg	SL267381

## Hygiene with Head Support



Small	900	190	570	620	995	290	440	850	205 kg	SL246181
Medium	1030	190	580	690	1010	290	440	980	205 kg	SL246281
Large	1100	190	630	800	1100	290	440	1050	300 kg	SL246381
Extra Large	1300	200	700	880	1230	300	450	1200	300 kg	SL246481
Extra Extra Large	1450	210	750	980	1370	300	450	1300	300 kg	SL246581

It is very important to use the correct sized sling and ensure it is properly fitted before attempting to lift. This will ensure the patient feels safe, dignified and comfortable, and allows the carer to be confident of an effective and safe procedure. Please refer to the ProSling Size Selection Chart for guidance on sizing recommendations. This chart is published in the ProSling Range Guide and is available to download from [www.novis.com.au](http://www.novis.com.au)

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.