

FOR PROFESSIONALS

Cochlear™ Baha® System

Product catalogue



Hear now. And always





Contents

-  1. SOUND PROCESSORS
-  2. ACCESSORIES
-  3. WIRELESS ACCESSORIES
-  4. SURGICAL PRODUCTS
-  5. NON-SURGICAL PRODUCTS
-  6. TRAINING AND DEMO PRODUCTS
-  7. FITTING SOFTWARE
-  8. APPS

Sound processors

BAHA® 5 SOUND PROCESSOR

Included in the order:

- Sound processor
- Abutment cover
- CS44 programming port covers
- Battery multi-pack
- Cleaning kit
- Safety line
- Test rod
- Battery door tool
- Document pack



reddot award 2015
winner

When ordering, specify your language requirements for the document pack. It contains important information such as the User Manual and the Warranty Card.



95201 Blonde



95202 Black



95203 Silver



95204 Brown



95205 Copper

Not all products are available on all markets.

Sound processors

BAHA 5 POWER SOUND PROCESSOR

Included in the order:

- Sound processor
- Abutment cover
- CS44 programming port covers
- Battery multi-pack
- Cleaning kit
- Safety line
- Test rod
- Tamper proof battery door tool
- Document pack

When ordering, specify your language requirements for the document pack. It contains important information such as the User Manual and the Warranty Card.



95470 Blonde



95471 Black



95472 Silver



95473 Brown



95474 Copper

Not all products are available on all markets.

Sound processors

BAHA 5 SUPERPOWER SOUND PROCESSOR

When ordering, specify your language requirements for the document pack. It contains important information such as the User Manual and the Warranty Card.



* Processing Unit

Not all products are available on all markets.

Sound processors

BAHA 5 SUPERPOWER SOUND PROCESSOR

Sound processing unit*

* Processing Unit



96001 Carbon



96002 Smoke



96003 Mocha



96004 Maize

Actuator unit



95461 Carbon



95462 Smoke



95463 Mocha



95464 Maize

Not all products are available on all markets.

Accessories

BAHA 5 SOUND PROCESSOR

Battery doors



P1129017 Baha 5 Sound Processor Tamper-proof battery door (including 3 battery doors, tool and sticker sheet)



P740869 Baha 5 Sound Processor Battery door, 3 pcs (including 3 battery doors and tool)



95330 Baha 5 Sound Processor Battery door (including 10 battery doors and tool)

Case

Not for US market



95224 Baha 5 Sound Processor Inner Case

Battery

Not for US market

95223 Battery, size 312, mercury free. (1 box of 10 cards, 1 card includes 6 batteries)

Not all products are available on all markets.

Accessories

BAHA 5 POWER SOUND PROCESSOR

Battery door



95478 Baha 5 Power Battery door
(including 3 battery doors, tool and sticker sheet)

Case

Not for US market.



95980 Baha 5 Power
Inner case

Battery

Not for US market.

94738 Battery, size 675, mercury free.
(1 box of 10 cards, 1 card includes 6 batteries)

Not all products are available on all markets.

Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

Rechargeable batteries



Standard



Compact



Carbon

Z285983 Standard

Z285988 Compact



Smoke

Z285984 Standard

Z285989 Compact



Mocha

Z285986 Standard

Z285991 Compact



Maize

Z285985 Standard

Z285990 Compact

Charging



Z341685

Battery Charger



Z342011

Charging Kit

Not all products are available on all markets.


Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

Cables

These are the recommended cables for Baha 5 SuperPower:

- 6cm, 8cm, and 11cm cables are for Same ear fittings
- 28cm cables are for Other ear fittings
- 50cm cables are for Bodyworn fittings



Carbon

Z299503 6 cm
Z299507 8 cm
Z299515 11 cm
Z299519 28 cm
95938 50 cm



Smoke

Z299504 6 cm
Z299508 8 cm
Z299516 11 cm
Z299521 28 cm
95939 50 cm



Mocha

Z299506 6 cm
Z299513 8 cm
Z299518 11 cm
Z299525 28 cm
95937 50 cm



Maize

Z299505 6 cm
Z299512 8 cm
Z299517 11 cm
Z299522 28 cm
95936 50 cm

Not all products are available on all markets.

Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

Earhook

For Behind-the-ear (BTE).



Z285994
Small, 3 pcs
Fitting tool and 3 fitting pins included (tamper resistant)



Z285995
Medium, 3 pcs



Z285996
Large, 3 pcs

Snugfit

For Behind-the-ear (BTE).



Small



Medium



Large



Carbon
Z285998 Small
Z285999 Medium
Z286000 Large



Smoke
Z286001 Small
Z299502 Medium
Z299510 Large



Mocha
Z299520 Small
Z299524 Medium
Z299523 Large



Maize
Z299511 Small
Z299514 Medium
Z299509 Large

Not all products are available on all markets.

Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

LiteWear

For Clip-on (Body-worn).



Z327614 LiteWear Fixing Aids
Includes: Alligator clip, safety pin, hook & loop.

Z327613 LiteWear Case

Earhook+

For Under-the-ear (UTE).



Medium

Large



Carbon
P771900 Medium
P771904 Large



Smoke
P771901 Medium
P771905 Large



Mocha
P771902 Medium
P771906 Large



Maize
P771903 Medium
P771907 Large

Not all products are available on all markets.

Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

Processor covers



White, 2 pcs
Z286007 Standard
Z319157 Compact



Blue, 2 pcs
Z286008 Standard
Z319158 Compact



Red Swirls, 2 pcs
Z319151 Standard
Z319163 Compact



Smoke, 2 pcs
Z319148 Standard
Z319160 Compact



Green, 2 pcs
Z286004 Standard
Z319154 Compact



Orange, 2 pcs
Z286005 Standard
Z319155 Compact



Seashells, 2 pcs
Z319149 Standard
Z319161 Compact



Blue Swirls, 2 pcs
Z319152 Standard
Z319164 Compact



Circles, 2 pcs
Z319150 Standard
Z319162 Compact



Pink, 2 pcs
Z319147 Standard
Z319159 Compact



Lavender, 2 pcs
Z286006 Standard
Z319156 Compact



Skulls, 2 pcs
Z319153 Standard
Z319165 Compact

Not all products are available on all markets.

Accessories

BAHA 5 SUPERPOWER SOUND PROCESSOR

Other accessories



Z339244 Adhesive labels, 15 pcs



Z300304 Mic Lock



Microphone Protector
Applicator with Removal
Tool

Z322658 2 pcs
Z354058 4 pcs



Z325644 Dry aid kits,
Europe
Z60761 Dry aid kits,
Australia



Z60764 Dry Brik by Dry
& Store®, 3 pcs

Not valid in US market.



Z342010 Storage case

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Extra battery doors: Baha 4, BP100



92079 Champagne blonde



92084 Soft black



92083 Slate grey



92082 Chestnut brown



92081 Glacier white



92080 Piano black



92572 Ruby red



92574 Iris purple



92573 Ocean blue



92575 Apple green

Battery: Baha 4, BP100

94330 Battery, size 13, mercury free.
(1 box of 100 cards, 1 card includes 6 batteries.)

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Extra battery doors: BP110 Power



92832 Champagne blonde



92833 Soft black



92834 Slate grey



92835 Chestnut brown



92836 Glacier white



92837 Piano black



93089 Ruby red



93090 Iris purple



93091 Ocean blue



93092 Apple green

Battery: BP110 Power

94331 Battery, size 675, mercury free.
(1 box of 100 cards, 1 card includes 6 batteries.)

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Tamper-proof battery doors: Baha 4, BP100

Included in the order:

- Tamper-proof battery doors
- Lock/unlock tool
- Sticker sheet



92077 Champagne blonde



92078 Soft black



92570 Peony pink



92571 Sky blue

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Tamper-proof battery doors: BP110 Power

Included in the order:

- Tamper-proof battery doors
- Lock/unlock tool
- Sticker sheet



92830 Champagne blonde



92831 Soft black



93088 Peony pink



93087 Sky blue

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Audio accessories

CAUTION: Always use the audio adapter to connect mains supplied electrical equipment to the sound processor.

For Cordelle II

- 90623** Audio cord
Beige, 75 cm
- 90625** Audio cord
Beige, 150 cm
- 90627** Audio cord
Black, 75 cm
- 90629** Audio cord
Black, 150 cm

For Cordelle II

- 90067** Audio adapter

Battery and Battery charger



- 90328** Battery charger 2x9V, EUR-plug
- 90332** Battery charger 9V, UK-plug
- HCB 663** Battery charger 9V, US-plug.
Used for recharging the 9V Cordelle II batteries.

- 95747** Battery size 9V.
For Cordelle II.

Actuator unit Baha Cordelle II



- 90015** Beige
- 90016** Black
- 90017** Grey

Not all products are available on all markets.

Accessories

LEGACY SOUND PROCESSORS

Audio accessories

CAUTION: Always use the audio adapter to connect mains supplied electrical equipment to the sound processor.

For BP110 Power and Baha 4

91868 Cochlear Baha Audio adapter long cable, 5 m.

For BP110 Power and Baha 4

91180 Cochlear Baha Audio adapter short cable, 0.7 m.

For BP110 Power and Baha 4

91181 Cochlear Baha Telecoil.

Not all products are available on all markets.

Accessories

BAHA ATTRACT SYSTEM

SP Magnet and Baha SoftWear™ Pads*

The SP Magnet includes Baha SoftWear™ Pads* (6 pcs).

*SP Magnet soft pads.



95771 SP Magnet 1



95772 SP Magnet 2



95773 SP Magnet 3



95774 SP Magnet 4



95775 SP Magnet 5



95776 SP Magnet 6



94175 SP Magnet soft pads, 24 pcs
Can be ordered separately.

Colour covers

The Colour covers includes a centering pin for easy application.



95190 Blonde, 3 pcs



95191 Brown, 3 pcs



95192 Silver, 3 pcs



95193 Copper, 3 pcs



96005 Maize, 3 pcs

Case



94125 End user case

Not all products are available on all markets.

Accessories

OTHER ACCESSORIES

Abutment covers



90744 Blonde



90745 Black



90746 Silver grey



90747 Brown

Miscellaneous



90711 Safety line



93379 Stickers
(3 sheets per pack)



95458 Baha Sound
processor Cleaning kit



92993 Accessory box



P1316179 Cochlear
Magnetic Battery
Removal Tool



92657 Screwdriver for
Tamperproof battery
doors

Not all products are available on all markets.

Wireless accessories



Cochlear Wireless Mini Microphone 2
P770841 GB
P770842 EU
P770844 US



Cochlear Wireless Phone Clip
94770 EU
94771 US
94772 GB
94773 AUS



Cochlear Baha Remote Control 2
94790 EU
94791 US
94792 GB
94793 AUS



Cochlear Wireless Mini Microphone 2+
P770845 GB
P770846 EU
P770847 AUS
P770848 US



Cochlear Wireless TV Streamer
94760 EU
94761 US
94762 GB
94763 AUS



94991 Cochlear Wireless demo kit
Delivered empty.

Not all products are available on all markets.

Surgical products

IMPLANTS

Cochlear Baha BI300 Implants and cover screw*

Used for Baha Connect two-stage surgery and Baha Attract surgery.

NOTE: The cover screw conical is only compatible with the BI300 Implants.



92128 BI300
Implant 3 mm



92129 BI300
Implant 4 mm



92136 Cover screw
conical

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

BAHA CONNECT SYSTEM

Cochlear Baha Implants with DermaLock™ Abutments*

Used for Baha FAST surgery.



93329 BIA400
Implant 4 mm
Abutment 6 mm



93330 BIA400
Implant 4 mm
Abutment 8 mm



93331 BIA400
Implant 4 mm
Abutment 10 mm



93332 BIA400
Implant 4 mm
Abutment 12 mm



93338 BIA400
Implant 4 mm
Abutment 14 mm

Cochlear Baha DermaLock Abutments*

Used for Baha two-stage surgery.

NOTE: The DermaLock Abutments (BA400) are only compatible with the BI300 Implants.



93333 BA400
Abutment 6 mm



93334 BA400
Abutment 8 mm



93335 BA400
Abutment 10 mm



93336 BA400
Abutment 12 mm



93337 BA400
Abutment 14 mm

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

BAHA CONNECT SYSTEM

Cochlear Baha BIA300 Implants with Abutments*

Used for Baha FAST surgery.



92126 BIA300
Implant 3 mm
Abutment 6 mm



92127 BIA300
Implant 4 mm
Abutment 6 mm



92346 BIA300
Implant 4 mm
Abutment 9 mm

Cochlear Baha BA300 Abutments*

The BA300 Series abutments are only compatible with the BIA300 implants.



92130 BA300
Abutment 6 mm



92131 BA300
Abutment 9 mm



93357 BA300
Abutment 12 mm

Cochlear Baha BA210 Abutments*

NOTE: The BA210 Series abutments are compatible with the previous Baha implants with hexagonal connection.



92132 BA210
Abutment 5.5 mm for
flange fixture
For hexagonal
connection.



92133 BA210
Abutment 8.5 mm for
flange fixture
For hexagonal
connection.

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

BAHA CONNECT SYSTEM

Cochlear Baha BIA310 Implants with Abutments*

Used for Baha FAST surgery.



P1340888 BIA310
Implant 4 mm with
6 mm abutment



P1340889 BIA310
Implant 4 mm with
8 mm abutment



P1340890 BIA310
Implant 4 mm with
10 mm abutment



P1340891 BIA310
Implant 4 mm with
12 mm abutment



P1340893 BIA310
Implant 4 mm with
14 mm abutment

Cochlear Baha BA310 Abutments*

The BA310 Series abutments are only compatible with the BI300 implants.



P1340894 BA310
Abutment 6 mm



P1340895 BA310
Abutment 8 mm



P1340896 BA310
Abutment 10 mm



P1340897 BA310
Abutment 12 mm



P1340898 BA310
Abutment 14 mm

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

BAHA ATTRACT SYSTEM

Cochlear Baha BIM400 Implant Magnet*

NOTE: The BIM400 Implant Magnet is only compatible with the BI300 Implants.



93550 BIM400 Implant Magnet

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

DISPOSABLE INSTRUMENTS

Healing caps*



95083 Healing cap
with plug Ø 20 mm



95084 Healing cap
with plug Ø 30 mm

Biopsy punch*



94870 Biopsy punch Ø 5 mm (10 pcs)



90940 Biopsy punch 4mm

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

DISPOSABLE INSTRUMENTS

Surgical templates for Baha 5 SuperPower Sound Processor*

For US market.



For US market:

P806653

Indicator for Baha 5
SuperPower, Baha
Connect System



For US market:

P806693

Indicator for Baha 5
SuperPower, Baha
Attract System

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

DISPOSABLE INSTRUMENTS

Drills*

NOTE: The Conical guide drill has a removable spacer for drilling to different depths.



93363 Conical guide drill 3+4 mm



92141 Widening drill 4 mm with countersink



92140 Widening drill 3 mm with countersink

* All sterile products are for single use only.

Not all products are available on all markets.

Surgical products

REUSABLE INSTRUMENTS

Baha instrument organiser

NOTE: The Baha instrument organiser includes the Component tray (92145), but does not include any instruments.



92146 Baha instrument organiser



92145 Component tray

Baha Ruler



93339 Baha Ruler

Not all products are available on all markets.

Surgical products

REUSABLE INSTRUMENTS

Surgical templates for Baha 5 SuperPower Sound Processor

Not for US market.



95934 Indicator for
Baha 5 SuperPower,
Baha Connect System



96048 Indicator for
Baha 5 SuperPower,
Baha Attract System

Not all products are available on all markets.

Surgical products

REUSABLE INSTRUMENTS

92147 Baha instrument set

The set includes the instruments to the right.

The instruments can also be ordered individually.

For cleaning and sterilisation guidelines see the Baha Support Guide for Surgery.



91095 Indicator for Baha



91116 Drill indicator for WS-75 and Osscora



90469 Screwdriver Unigrip 95 mm



92143 Multi wrench with ISO adapter



90381 Machine screwdriver Unigrip 25 mm



92142 Implant inserter



HIA 009-0 Counter torque wrench



90944 Raspatorium



90943 Dissector



90474 Abutment inserter

Wrench adapters



92347 ISO adapter for Multi wrench.
Spare part for Multi wrench.
For compatibility with ISO instruments such as implant inserter.

Not all products are available on all markets.

Surgical products

REUSABLE INSTRUMENTS

Baha Attract surgical instruments

NOTE: The instruments are used for the Baha Attract System surgical procedure, with BI300 Implant and BIM400 Magnet.



93571 Indicator for Baha Attract



94071 Implant Magnet template



95180 Bone bed indicator



95070 Soft tissue gauge 6 mm

Not all products are available on all markets.

Surgical products

OSSCORA ACCESSORIES AND SPARE PARTS

Accessories



91061 Irrigation tubing set 3.8 m, 6 pcs, each tubing is single use only



91102 Sterilisation cassette



91051 Osscora Surgical set 230V (Europe)
91052 Osscora Surgical set 230V (Asia Pacific)
91053 Osscora Surgical set 115V

Spare parts



91050 Osscora handpiece 20:1



91054 Osscora motor with 3.5 m cable



91891 Spray clip, 3 pcs (for Osscora handpiece)



04032600 Motor protective cover



92106 Clips for Osscora motor cable, 5 pcs

Not all products are available on all markets.

Non-surgical products

BAHA SOUNDARC

The Baha SoundArc is delivered as unilateral and it is possible to make it bilateral by adding a second Connector disc (P807181).



Baha SoundArc

P1324401 Extra small

P807175 Small

P807176 Medium

P807177 Large

P807178 Extra large



P807181 Baha SoundArc
Connector disc



95809 Baha Softband
soft pad, 6 pcs



P1572994

Baha SoundArc Hair matching Colour Kit

The article includes one pair of tips for each colour.



P1573011

Baha SoundArc Fun Colour Kit

The article includes one pair of tips for each colour.

Not all products are available on all markets.

Non-surgical products

BAHA SOFTBAND

The Softband is available in both unilateral and bilateral models. It comes in one adjustable size to fit all.



Unilateral



Bilateral



Pink
95750 Unilateral
95760 Bilateral



Light blue
95751 Unilateral
95761 Bilateral



Dark blue
95752 Unilateral
95762 Bilateral



Brown
95753 Unilateral
95763 Bilateral



Black
95754 Unilateral
95764 Bilateral



Blonde
95755 Unilateral
95765 Bilateral



Red
95756 Unilateral
95766 Bilateral



Jungle
95757 Unilateral
95767 Bilateral



95809 Baha Softband
soft pad, 6 pcs

Not all products are available on all markets.

Training and demo products

Skull Models



91153 ½ Skull model with soft tissue



90643 ¼ Skull model, plastic temporal bone



90644 ½ Skull model, plastic temporal bone



P1571591 ½ Surgical skull model with soft tissue

Ear Model



P1607869
Cochlear Ear Model

Not all products are available on all markets.

Training and demo products

BAHA 5 SOUND PROCESSOR

Demo sound processor



95207 Baha 5 Demo

Dummy sound processor*

Dummy processor, non functioning.



95820 Baha 5 Dummy

Demonstration equipment*

NOTE: For information about the Baha 5 SuperPower Demo kit please contact your local Cochlear representative.



95208 Baha 5 and Baha 5 Power Colour Guide



95206 Baha 5 Systems Presentation Kit (includes Baha 5 Sound Processor Dummy, SP Magnet Dummy and Baha 5 Colour guide)

* Not for clinical use.

Training and demo products

BAHA 5 POWER SOUND PROCESSOR

Demo sound processor



95475 Baha 5 Power
Demo

Dummy sound processor*

Non functioning.



95476 Baha 5 Power
Dummy

* Not for clinical use.

Not all products are available on all markets.

Training and demo products

BAHA 5 SUPERPOWER SOUND PROCESSOR

Demo sound processor



96000 Baha 5
SuperPower Processing
Unit Demo

Dummy sound processor*

Non functioning.



95465 Baha 5
SuperPower Dummy

Demonstration equipment

NOTE: For information about the Baha 5 SuperPower Demo kit please contact your local Cochlear representative.



96051 Baha 5
SuperPower Colour
Guide

* Not for clinical use.

Not all products are available on all markets.

Training and demo products

Dummy Implants and Abutments*



92150 Dummy BI300
Implant 4 mm*



92149 Dummy BIA300
Implant 4 mm with
Abutment 6 mm*



93362 Dummy BIA400
Implant 4 mm with
Abutment 10 mm*



93426 King Size
Dummy Implant 4 mm
with Abutment 10 mm*



P1353482
Dummy BA310
Abutment 10mm



P1353483
Dummy BIA310 4mm
Implant with 10mm
abutment



P1416164
King Size Dummy
Implant 4 mm w BA310
Abutment 10 mm

Test Rod

NOTE: For information about the Baha 5 SuperPower Demo kit please contact your local Cochlear representative.



92612 Test Rod

Demo Baha Connect System*



92151 Snap trainer
BIA300

* Not for clinical use.

Not all products are available on all markets.

Training and demo products

BAHA ATTRACT SYSTEM

Dummy Magnets for Baha Attract System*

NOTE: The Dummy BIM400 Implant Magnet does not include an implant.



93580 Dummy BIM400 Implant Magnet*



95821 SP Magnet dummy*

Demo box Baha Attract System

Delivered empty.

NOTE: The SP Magnets, the accessory box (92933) and any accessories are ordered separately.



95152 Baha Attract SP Magnet selection kit box

* Not for clinical use.

Not all products are available on all markets.

Fitting software

BAHA FITTING SOFTWARE 5.4

Baha Fitting Software – Cochlear Software Distribution System (CSDS)

Please contact your local Cochlear representative to receive the link to the CSDS to download the Baha Fitting Software.

Available via download:

- Baha Fitting Software 5.4, 5.4 KEY, or 5.4 Basic
- BFS Installation Instruction
- BFS User Guide

Baha Fitting Software – CD

Available on CD:

- Baha Fitting Software 5.4 or 5.4 KEY
- BFS Installation Instruction
- BFS User Guide

P803254 Baha Fitting Software
(BFS 5.4, software with full access to all products and functions)

95430 Baha Fitting Software
(BFS 5.4 KEY, software with restricted access requiring a product key to unlock new products or features. Please find the sales numbers for the product keys on the next page.)

Programming cables

These articles can be ordered as spare parts when fitting a:

- Baha 5 Sound Processor
- Baha 5 Power Sound Processor
- Baha 5 SuperPower Sound Processor*

* *Baha 5 SuperPower Fitting adaptor required*



95948 Baha 5
SuperPower Fitting
adaptor
(9 cm, 3.5 in long)



95386 Baha
Programming cable
CS44 short blue
(50 cm, 19.7 in long)



95387 Baha
Programming cable
CS44 short red
(50 cm, 19.7 in long)



95388 Baha
Programming cable
CS44 long blue
(200 cm, 78.7 in long)



95389 Baha
Programming cable
CS44 long red
(200 cm, 78.7 in long)

Not all products are available on all markets.

Fitting software

BAHA FITTING SOFTWARE 5.4

Additional

Printed manuals (order separately)

Please contact your local Cochlear representative to receive a printed version of the material you require, or download the manuals at www.cochlear.com/manuals

Promotional material

Baha Fitting Software – Additional promotional material, including Desk Reference. (Contact your local representative.)

Product keys

The product keys can be used if a restricted software version has been installed previously. The specific product key will only provide a specific type of access.

NOTE: As soon as a new product key is being entered, previous accessed content will be overwritten.

Baha Generation 5 product keys:

P1325948 Product key-access to SoundArc

P1325955 Product key-access all Baha 5 generation sound processors, excluding BA300 connection type

P1325949 Product key-access all Baha 5 generation sound processors, excluding BA300 and SoundArc connection types

P1325950 Product key-access all Baha 5 generation sound processors excluding SuperPower, excluding BA300 connection type

P1325951 Product key-all Baha 5 generation sound processors and connection types excluding SuperPower, BA300 and SoundArc

P1325952 Product key-all Baha 5 generation sound processors and connection types excluding SuperPower, DermaLock

P1325953 Product key-all Baha 5 generation sound processors and connection types excluding SuperPower, DermaLock, SoundArc

P1325954 Product key-access to BA300 and SoundArc

Due to local regulatory requirements not all content is available in all markets.

Please contact your Cochlear representative to make sure you order the appropriate product key.

Not all products are available on all markets.

Fitting software

BAHA FITTING SOFTWARE 4.0 SR3

Baha Fitting Software 4.0 SR3 – Cochlear Software Distribution System (CSDS)

Please contact your local Cochlear representative to receive the link to the CSDS to download the Baha Fitting Software 4.0 SR3.

Available via download:

- Baha Fitting Software 4.0 SR3 and 4.0 SR3 Basic
- BFS Installation Instruction
- BFS User Guide

94227 Baha Fitting Software – CSDS (Download)

(Software with full access to all products and functions. Please find the sales numbers for the product keys on the next page.)

Baha Fitting Software 4.0 SR3 – CD

Available on CD:

- Baha Fitting Software 4.0 SR3
- BFS Installation Instruction
- BFS User Guide

94223 Baha Fitting Software 4.0 SR3

(Software with full access to all products and functions. Please find the sales numbers for the product keys on the next page.)

Not all products are available on all markets.

Fitting software

BAHA FITTING SOFTWARE 4.0 SR3

Additional

Printed manuals (order separately)

Please contact your local Cochlear representative to receive a printed version of the material you require, or download the manuals at www.cochlear.com/manuals.

Promotional material

Baha Fitting Software – Additional promotional material, including Desk Reference. (Contact your local representative.)

Product keys

The articles can be used if a restricted software version has been installed previously. The specific product key will only provide a specific type of access.

For use with BFS 4.0 SR3 Basic.

NOTE: As soon as a new product key is being entered, previous accessed content will be overwritten.

Baha Generation 3 and 4 product keys:

P1325958 Product key-access to generation 3 and 4 sound processors, Attract and wireless

P1325956 Product key-access all generation 3 sound processors and Attract

P1325957 Product key-access all generation 3 sound processors

Due to local regulatory requirements not all content is available in all markets.

Please contact your Cochlear representative to make sure you order the appropriate product key.

Not all products are available on all markets.

Apps

Baha 5 Smart App

The Baha 5 Smart App makes it possible to control and customise the Baha 5 sound processor directly from recipients' iPhone or Android smartphone. The app also provides functionality for personal preferences and customisation.

The Cochlear Baha 5 Smart App is available on App Store and Google Play and may not be available in all countries. For compatibility information visit www.cochlear.com/compatibility.

Available to download from Cochlear.com:

- Baha 5 Smart App Setup Guide

Compatible with:

- Baha 5 Sound Processor
- Baha 5 Power Sound Processor
- Baha 5 SuperPower Sound Processor



Baha 5 Smart



Not all products are available on all markets.

Apps

Baha Control App

The Baha Control App makes it possible to control the Baha sound processor from recipients' iPhone or Android smartphone via the Cochlear Wireless Phone Clip.

The Cochlear Baha 5 Smart App is available on App Store and Google Play and may not be available in all countries. For compatibility information visit www.cochlear.com/compatibility.

Available to download from Cochlear.com:

- Baha Control App Setup Guide

Compatible with:

- Baha 5 Sound Processor
- Baha 5 Power Sound Processor
- Baha 5 SuperPower Sound Processor
- Baha 4 Sound Processor

Required device:

- Cochlear Wireless Phone Clip



Baha Control



Not all products are available on all markets.

Apps

Baha Support App

The Baha Support App is created to provide information and assistance for Baha recipients. It contains instructions for use and short video tutorials to get the most out of the Baha System.

The Cochlear Baha 5 Smart App is available on App Store and Google Play and may not be available in all countries. For compatibility information visit www.cochlear.com/compatibility.

Contains support for:

- Baha 5 Sound Processor
- Baha 5 Power Sound Processor
- Baha 5 SuperPower Sound Processor
- Baha 4 Sound Processor
- Baha BP110 Power Sound Processor
- Baha BP100 Sound Processor



Baha Support




Not all products are available on all markets.

Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to helping people with moderate to profound hearing loss experience a life full of hearing. We have provided more than 550,000 implantable devices, helping people of all ages to hear and connect with life's opportunities.

We aim to give people the best lifelong hearing experience and access to innovative future technologies. We have the industry's best clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

 **Cochlear Bone Anchored Solutions AB**, Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden
Tel: +46 31 792 44 00, Fax: +46 31 792 46 95

Regional Offices

Cochlear Ltd, (ABN 96 002 618 073), 1 University Avenue, Macquarie University, NSW 2109 Australia
Tel: +61 2 9428 6555, Fax: +61 2 9428 6352

Cochlear Americas, 13059 E. Peakview Avenue, Centennial, CO 80111, USA
Tel: +1 303 790 9010, Fax: +1 303 792 9025

Local Offices

Cochlear Europe Ltd, 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey, KT15 2HJ, United Kingdom
Tel: +44 1932 26 3400, Fax: +44 1932 26 3426

Cochlear Deutschland GmbH & Co. KG, Karl-Wiechert-Allee 76A, 30625 Hannover, Germany
Tel: +49 511 542 770, Fax: +49 511 542 7770

Cochlear Benelux NV, Schaliënhoedreef 20 I, 2800 Mechelen, Belgium
Tel: +32 15 79 55 77, Fax: +32 15 79 55 70

Cochlear France S.A.S., 135 route de Saint Simon, 31035 Toulouse, France
Tel: +33 5 34 63 85 85 (international), Tel: 0805 200 016 (national), Fax: +33 5 34 63 85 80

Cochlear Italia SRL, Via Larga 33, 40138 Bologna, Italy
Tel: +39 051 601 53 11, Fax: +39 051 392 062

Cochlear Tıbbi Cihazlar ve Sağlık Hizmetleri Ltd. Sti., Cubuklu Mah. Bogazici Cad. Bogazici Plaza, No: 6/1 Kavacik, TR-34805 Beykoz-Istanbul, Turkey
Tel: +90 216 538 59 00, Fax: +90 216 538 59 19

Cochlear Nordic AB, Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden
Tel: +46 31 335 14 61, Fax: +46 31 335 14 60

Cochlear Canada Inc, 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 416 972 5082, Fax: +1 416 972 5083

Nihon Cochlear Co Ltd, Ochanomizu-Motomachi Bldg 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan
Tel: +81 3 3817 0241, Fax: +81 3 3817 0245

Cochlear AG, EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland
Tel: +41 61 205 8204, Fax: +41 61 205 8205

Cochlear Latinoamerica, S. A., International Business Park Building 3835, Office 403 Panama Pacifico, Panama
Tel: +507 830 6220, Fax: +507 830 6218

Cochlear Limited (Singapore Branch), 238A Thomson Road #25-06, Novena Square Office Tower A, Singapore 307684, Singapore
Phone: +65 65533814, Fax: +65 64514105

Cochlear Medical Device (Beijing) Co Ltd, Unit 2608-2617, 26th Floor, No.9 Building, No.91 Jianguo Road, Chaoyang District, Beijing 100022, PR China
Tel: +86 10 5909 7800, Fax: +86 10 5909 7900

Cochlear (HK) Ltd, Room 1404-1406, 14/F, Leighton Centre, 77 Leighton Road, Causeway Bay, Hong Kong
Tel: +852 2530 5773, Fax: +852 2530 5183

Cochlear Korea Ltd, 1st floor, Cheongwon Building 33, Teheran-ro 8 gil, Gangnam-gu, Seoul, Korea
Tel: +82 2 533 4450, Fax: +82 2 533 8408

Cochlear Medical Device Company India PVT Ltd, Platina Bldg, Ground Floor, Plot No. C 59, G Block, BKC, Bandra East, Mumbai 400051 India
Tel: +91 22 6112 1111, Fax: +91 22 61121100

Cochlear Colombia, Avenida Carrera 9 #115-06 Of. 1201 Edificio Tierra Firme, Bogota D.C., Colombia
Tel: +57 315 339 7169 / +57 315 332 5483

Cochlear México, S.S.A. de C.V. Av. Tamaulipas 150 Torre A piso 9, Col. Hipódromo Condesa, 06170 Cuauhtémoc, Ciudad de México, México
Tel: +52 0155 5256 2199

www.cochlear.com

This content is meant for professional use. If you are a consumer, please seek advice from your medical practitioner or health professional about treatments for hearing loss. They will be able to advise you on a suitable solution for your hearing loss condition. All products should be used only as directed by your medical practitioner or health professional. Not all products are available in all countries. Please contact your local Cochlear representative.

Cochlear, Baha, 科利耳, コクレア, 코클리어, Hear now. And always, SmartSound, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB or Cochlear Limited (unless otherwise noted). Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android is a trademark of Google LLC. Google Play and the Google Play logo are trademarks of Google Inc.

© Cochlear Bone Anchored Solutions AB 2019. All rights reserved. DEC19.

D1309440-V7

