

# Bovie®

OPERATING WITH PRECISION



## QUALITY SURGICAL PRODUCTS, WITH PROCEDURE-FOCUSED DESIGN

Since its establishment in the late 1970s, Bovie Medical has been striving to continue with the pioneering work of Dr. Bovie, credited as the inventor of modern electrosurgery, by continuing to develop and manufacture advanced electrosurgical products. Manufactured in the USA, Bovie Medical is recognised for their quality design and manufacturing as well as excellent customer service. It has become a leading manufacturer that sets the standard for electrosurgical generators and accessories worldwide.

## BOVIE HIGH FREQUENCY DESICCATORS - IT'S YOUR CHOICE




### Advantages of choosing Bovie;

Three, low-end electrosurgical devices available  
Choice of modes for varying applications  
Low-cost models for everyday procedures  
Standard accessories across the range  
Huge choice of dermal tips for a wide range of procedures  
Dermal tips are compatible with (and often cheaper than),  
the Conmed Hyfrecator® product  
Made in the USA with a 4-year warranty



ELECTROSURGERY - BETTER BY BOVIE!

# BOVIE'S SUITE OF HIGH FREQUENCY DESICCATORS FOR OFFICE USE

	DERM 102	DERM 942	BANTAM PRO A952
			
<b>Overview</b>	The Derm 102 for Family Practice is an economical, high frequency desiccator, designed to give practitioners the ability to perform minor skin procedures in-office.	Derm 942 - designed and engineered to be the most reliable and durable 40 watt desiccator available today for a wider range of applications.	The Bovie Bantam PRO, A952 is a High Frequency Desiccator with all the capabilities of the Derm 942, plus Cutting ability.
<b>Max. Power</b>	10 watts	40 watts	50 watts (40W Bipolar)
<b>Adjustment</b>	0.1W	0.1W (up to 10W) then 1W	0.1W (up to 10W) then 1W
<b>Monopolar</b>	Yes	Yes	Yes
<b>Bipolar</b>	Yes	Yes	Yes
<b>Warranty</b>	4 years	4 years	4 years

Products and specifications are subject to change without notice. ©2023 Bovie Medical and DocStock.

## Common Office Procedures

DESICCATION	LESION	SUGGESTED POWER RANGE	LESION	SUGGESTED POWER RANGE	HEMOSTASIS		SUGGESTED POWER RANGE
	Angiomas (Cherry)	1.2 - 2.0 watts	Pyogenic Granuloma	5.0 - 10.0 watts		Ellipse excision	5.0 - 10.0 watts
Telangiectasias	1.5 - 2.0 watts	C and D	5.0 - 10.0 watts	Punch excision	5.0 - 10.0 watts		
Seborrheic Keratosis	3.0 - 4.0 watts	Skin Tags (Acrochordons)	1.0 - 3.0 watts	Shave excision	5.0 - 10.0 watts		
Warts	5.0 - 15.0 watts						

## Accessories

Bovie Medical provides a complete portfolio of accessory options for your delicate procedural needs. Below are examples of most commonly used accessories. For a more complete list, please visit our website - [www.DocStock.com.au](http://www.DocStock.com.au)



Sharp Disposable Tips - Sterile/Non Sterile  
BV-A805/BV-A804



Blunt Disposable Dermal Tips - Sterile/Non Sterile  
BV-A806/BV-A807



Bipolar Forceps  
BV-381-250A/BV-381-250B/  
BV-381-250C



Reusable Angled Fine Tip Tungsten Epilation Needle  
BV-A34T



Disposable Handpiece Sheaths  
BV-A910



3-button Power Control Pencil and Lead  
BV-A902



Bipolar Cable for DERM 102  
BV-A827V



Footswitch  
BV-A803



Desk Bracket  
BV-A813



Mobile Stand and Bracket  
BV-1580ARN



**docStock**  
.com.au

**DocStock Pty Ltd**

ABN 48 115 268 203

Unit 20, 7 Anella Ave, Castle Hill NSW 2154  
Phone: (02) 8850 7177 Fax: (02) 8850 7988  
Email: [info@docstock.com.au](mailto:info@docstock.com.au)  
[www.DocStock.com.au](http://www.DocStock.com.au)

