

# ACEEVAC PORTABLE SUCTION PUMP

## MAINS OPERATED

ACeevac continues the tradition of the Clements Emergency pump, providing the same dependable level of performance in a lighter, more compact form factor.

Designed for high vacuum performance, reliability and minimal maintenance, ACeevac is ideal for clinical and home use or as an alternative to wall suction.

- Maximum Vacuum: -80 kPa [ -600 mmHg]
- Maximum Flow : 21 litres/minute
- Direct vacuum regulation
- Compact and lightweight

Supplied ready to use with:

- 1L Reusable jar x1
- Silicone 8 x 14 mm suction tubing, 150 cm x1, 18 cm x2
- Hydrophobic / antibacterial filter x1
- Stepped conical tubing connector x1



# CLEMENTS

To order or for more information, please contact

**docStock**  
MEDICAL EQUIPMENT SPECIALISTS™

DocStock Pty Ltd  
ABN 48 15 268 203  
P. (02) 8850 7177  
E. info@docstock.com.au  
www.DocStock.com.au

# ACEEVAC PORTABLE SUCTION PUMP

SPECIFICATIONS	
Item	SUC81025
Description	High Vacuum High Flow Mains Operated Portable Suction Pump with 1 x 1L Reusable Collection Jar, Bacterial/Hydrophobic Filter and Silicone Suction Tubing Set
Power Rating	220 - 240V 50/60Hz 110 VA
Maximum Vacuum	-80 kPa, [-600 mmHg]
Maximum Flow Rate	21 L/min (through jar) 32 L/min (free air)
Gauge	40 mm diameter, Bourdon tube type, dual scale, graduated in kPa and mmHg
Gauge Range	0 to -100 kPa graduated at 5 kPa 0 to -760 mmHg graduated at 20 mmHg
Vacuum Regulator	Needle valve type
Pump Type	Piston - motor and pump combined
Motor Type	PSC AC motor
Collection Canister	1 Litre autoclavable, shatterproof polycarbonate jar Optional: 1 Litre disposable liner jar, 2 Litre disposable liner jar
Overflow Protection	Float valve mechanism / Hydrophobic filter
Power Indicator	Neon indicator shows power available
Dimensions	360L x 215W x 195H mm
Weight	3.6 kg
GMDN	63642
Device Class	Class IIA
Electrical Protection	Class II
Operation Mode	Continuous operation

ACCESSORIES / SPARES	
SUC81030001	Kit, 1L Reusable Collection Jar (Pack 1) (1L Reusable Jar, Lid, Float Assembly)
SUC81030005	Support Jar for FLOVAC Disposable 1L Liner
SUC81030009	Kit, Seal, Lid, 1L / 2L Reusable Jar (Pack 2)
SUC81030010	Kit, Lid Assembly for 1L / 2L Reusable Jar (Pack 2)
SUC81030012	Kit, Filter, Bacterial/Hydrophobic (Pack 10)
SUC81030013	Kit, Silicone Tubing Set, Joiner (Pack 1)
SUC81030121	Kit, 1L Disposable Liner (Pack 50)
SUC81030122	Kit, 2L Disposable Liner (Pack 50)
SUC81030132	Support Jar for FLOVAC Disposable 2L Liner
SUC84500033	Nipple, Yellow for FLOVAC Disposable Liner

All specifications are subject to change without notice. Please check the exact specification when ordering.

Sponsor:  
ICU Medical Australia Pty Ltd  
Unit U, 10 - 16 South St, NSW 2116, Australia  
Telephone: 1300 438 652 +61 2 9466 5300  
Email: [au\\_orders@icumed.com](mailto:au_orders@icumed.com)  
Web: [www.clements.net.au](http://www.clements.net.au)

 CA-MI srl  
Via Ugo La Malfa 13  
Frazione Pilastro  
43013 Langhirano (PR)  
Italy



**CLEMENTS**

Rev 3. © Copyright 2021 ICU Medical Australia Pty Ltd