



DL5.au

SUPER POLARIZED. With the world's first **variable polarization system** in a hybrid dermatoscope, DermLite® DL5 lets you uncover a whole new layer in handheld dermoscopy. Not only can you toggle between polarized and non-polarized modes, now you're able to visualize extremely superficial structures under parallel polarization and, as you scroll through the entire polarization range, optically explore a lesion's visual depths.

And on that note, you'll love DL5's optics, which are stunning up-close and outstanding at a distance: Its **32 mm lens system** with *true* 10x magnification reveals the supertiniest of details — better than ever before.

CROSS-POLARIZED



NON-POLARIZED







DermLite DL5 supports your clinical workflow unlike any other dermatoscope. It's the world's first hybrid dermatoscope with **365nm UV**, allowing you to highlight fluorescent structures — at full 10x magnification.

And, when screening your patients at a distance, DL5's **high-powered torch LED** provides phenomenal illumination to larger areas or hard-to-reach spots.

SUPER CLINICAL.

Need to measure a lesion? DL5 comes with an on-board, magnetically attached 100 mm ruler. Simply slide it off and its circular markings make it easier than ever to assess a lesion or choose the proper biopsy punch.





ALL IN.

DL5 comes with everything you need for visual and digital dermoscopy right out of the box.

> USB TO USC-C CABLE, 2M



SUPER Charged.

DL5 lets you work substantially longer on a single charge, thanks to a battery that's more than **twice the capacity of DL4**.

Between uses, keep it charged up on the included desktop charging base — which doubles as home base for your IceCap® supply. You'll even be able to charge a second DermLite or your phone using its front-facing 5V USB port. Or, use the supplied USB-C cable while on the go; regardless, DL5's 4-level charge indicator keeps you up to date.





- **POLARIZATION** Variable polarization (incl. polarized/non-polarized toggle)
 - **OPTICS** ø32 mm, 10x
- **ILLUMINATION** 16x white, 4x UV (365nm), 4x PigmentBoost, 1 torch LED (each mode is brightness-adjustable)
 - **POWER** 3000 mAh 3.7V lithium-ion, 4-level charge indicator, includes charging base with second 5V output
- **SMARTPHONE USE** Via included MCC[™] smartphone adapter
- **CONTACT PLATE** Removable via bayonet (incl. 10mm reticle)
- **INFECTION CONTROL** IceCap[®] disposable covers (5 included)
- INCLUDED ACCESSORIESMagnetic eyepiece, MCC smartphone adapter, 100 mm steelruler, leather belt pouch, protective holster, 5 IceCaps, 2 mUSB-A to USB-C charging cable, microfiber cleaning cloth

SIZE & MASS 182x58x36mm, 285g







©2022 by DermLite LLC. All rights reserved. For information about DermLite patents & trademarks, please visit www.dermlite.com/patents | DermLite LLC | 31521 Rancho Viejo Rd Ste 104 | San Juan Capistrano, CA 92675 USA Phone +1-949-481-6384 Fax +1-949-240-7492 info@dermlite.com | dermlite.com | Specifications subject to change without notice.

DL5.au