





- Non-touch measurement
- High precision infrared sensor
- Practical curved design
- High definition display



Alicn Non-touch Clinical Infrared Thermometer

A large, clear display and clean curved design makes this clinical forehead thermometer perfect for fast use. Take accurate readings within about 1 second. Features a single button operation, fever alert function, 32-unit memory and 2-yr warranty!



Specifications:

- Product Size : 40mm (L) x 44mm (W) x 164mm (H)
- Product Weight : 61g
- Memory : 32 groups
- Measuring mode : Forehead mode
- Measurement range : 32.0°C~42.2°C
- Accuracy : $\pm 0.2^{\circ}\text{C}$ (35.0°C~42.0°C); $\pm 0.3^{\circ}\text{C}$ (32.0°C~34.9°C / 42.1°C~42.2°C)
- Voltage source : 2 x 1.5 V Batteries: size AAA (included)
- Display : LED digital display
- Automatic power off : 60s
- Ordering code: AM-R1D2-LED

