

Unmatched Pulse Oximetry Performance and Value



WristOx₂ 3150 with BLE



WristOx₂ Sensor Line



WristOx₂ Accessories

WristOx₂® Model 3150 with Bluetooth Low Energy (BLE)

Wrist-worn Pulse Oximeter

Engineered with Nonin Medical's proven PureSAT® SpO₂ technology, the **WristOx₂ 3150 with BLE** delivers highly accurate readings in a wide range of patient populations and settings.

The **WristOx₂ 3150 with BLE** is ideal for applications including ambulatory monitoring, remote wireless monitoring and overnight studies—providing the versatility and reliability that clinicians expect from Nonin.

The most advanced wrist-worn pulse oximeter available, the **WristOx₂ 3150 with BLE** features improved performance and reporting capabilities.



For Proven Performance, Count on Nonin.

Experience Nonin Quality and Greater Affordability

- Small size wears like a watch with large display screen and enhanced sensor connector design
- Transfer data securely and wirelessly to the NoninConnect app*, available on the Apple Store or download directly to nVision via USB cable
- Excellent for many clinical applications such as the six-minute walk test and overnight sleep studies
- Dependable design and construction— with industry leading three year warranty



*The NoninConnect app isn't cleared as a medical device



Focus On What Matters Most

Clinical Efficiency



Memory Volume Indicator Mode provides verification of record data during oxygen qualification studies and overnight sleep studies



1,080 hours of non-volatile memory at four-second sampling rate



Nonin's WristOx₂ sensors have been designed exclusively for use with WristOx₂ pulse oximeters and includes a selection of three styles and sizes

Clinical Effectiveness



The WristOx₂ Model 3150 with BLE is an ideal choice for applications where small, portable oximetry is needed



Long lasting, 44 hour battery life



Available in two assemblies:
Single Unit: includes sensor + wristband
Starter Kit: includes sensor, three wristbands, USB download cable, and nVision data management

Nonin Medical, Inc.
13700 1st Avenue North
Plymouth, MN • 55441-5443 • USA

Tel: +1.763.553.9968 1.800.356.8874
Fax: +1.763.577.5521
Email: info@nonin.com

Nonin Medical B.V.
Doctor Paul Janssenweg 150
5026 RH Tilburg • Netherlands

Tel: +31 (0)13 – 45 87 130
Email: infointl@nonin.com



nonin.com



©2021 Nonin Medical, Inc. All trademarks are the property of Nonin Medical Inc. unless otherwise noted.
M-18-037-03

CE 0123



UL 60601-1
300M

Rx Only