

Patient Eye Protection


ANESTHESIA EYE PATCH

Product Reference # TT55

FUNCTION:

Protecting the eye during anesthesia by keeping the eye moist, preventing corneal abrasions and keeping particles out of the eyes.

<p>FEATURES:</p> <ul style="list-style-type: none"> • Secure Lock side Tab • Exclusive Vertical Removal Tab • Theia Patent Pending Comfort Zone • Pro4 Adhesive Ring • Hypoallergenic • Latex Free 	<p>STANDARDS *:</p> <ul style="list-style-type: none"> • Registered FDA Class 1 Medical Device • European Class 1 Medical Device • TT55 basic UDI-DI: 0710859306030
	<p>PACKAGING:</p> <ul style="list-style-type: none"> • 50 pair per box • Inner Box Dimensions: 2.75 in x 3.38 in x 3.25 in, 7 cm x 8.6 cm x 8.3 cm • Master Box Dimensions 18in x1 4in x1 0in, 45.7 cm x 35.6 cm x 25.4cm • 50 boxes per master carton, 50 Pair Per Box • Master Carton – 11 lbs, 5 kg • Box 50 pair - 3.5 oz, 100 gm



0 710859 306030


REF TT55

MD **CE**

! 10°C/50°C

Latex **Single Use**

THEIA EYE BLOCK
 3 Canal St.
 Cumberland RI 02865
 USA
 Tel (888) 626 3319
 www.theiaeye.com



THEIA EYE BLOCK

ANESTHESIA EYE SHIELD - 50PR/PK

REF TT55 **43877-1-1** **LOT**

THE-ANESTHESIA-G1 Rev. C

QTY - 50PK

50 PK

Jan/2023

Made in the U.S.A.