


## Patient Eye Protection **IPL EYE SHIELD** Product Reference # TT73

**FUNCTION:**

Provides protection from wavelengths emitted from IPL devices. The IPL eye shields block out the wavelengths, energy and extremely bright lights generated from IPL devices. Our unique layers actually reduce the energy by a factor of 10,000 while also shielding the bright lights from patients' eyes.

<p><b>FEATURES:</b></p> <ul style="list-style-type: none"> <li>Provides protection from IPL and LED medical devices</li> <li>Exclusive Removal Tab</li> <li>Theia Patent Pending Comfort Zone</li> <li>Pro4™ Adhesive Ring</li> <li>Hypoallergenic</li> <li>Latex Free</li> </ul>	<p><b>STANDARDS *:</b></p> <ul style="list-style-type: none"> <li>European Class 1 Medical Device</li> <li>Meets IPL ISO std 12609-1:2013-07</li> <li>Registered FDA Class 1 Medical Device</li> <li>Basic UDI-DI: 0710859306023</li> </ul>
	<p><b>PACKAGING:</b></p> <ul style="list-style-type: none"> <li>50 pair per box</li> <li>Inner Box Dimensions: 2.75 in x 3.38 in x 3.25 in, 7 cm x 8.6 cm x 8.3 cm</li> <li>Master Box Dimensions 18in x1 4in x1 0in, 45.7 cm x 35.6 cm x 25.4cm</li> <li>50 boxes per master carton, 50 pair per box</li> <li>Master Carton - 11lbs, 5kg</li> <li>Box 50 pair - 3.5oz, 99.8gm</li> </ul>




**REF TT73**



Single Use

**THEIA EYE BLOCK**  
 3 Canal St.  
 Cumberland RI 02865  
 USA  
 Tel (888) 626 3319  
 www.theiaeye.com



**THEIA EYE BLOCK**

**IPL-EYE SHIELD - 50PR/PK**

**REF TT73** **43877-1-1** **LOT**

**THE-IPL-G2 Rev. C**

**QTY - 50PK**

Jan/2023 **Made in the U.S.A.**