


Patient Eye Protection LASER EYE SHIELD Product Reference # TT63

FUNCTION:

Provides protection from all lasers operating between 315 nm - 11,000 nm with an optical density of 7. These shields also disperse the heat typically emitted from infrared lasers.

<p><u>FEATURES:</u></p> <ul style="list-style-type: none"> Exclusive Two-Tab design for easy placement and removal Theia Patent Pending Comfort Zone Pro4 Adhesive Ring Hypoallergenic Latex Free 	<p><u>STANDARDS *:</u></p> <ul style="list-style-type: none"> Meets damage threshold requirements of EN207 and ANSI Z 136.1 and .7 European Class 1 Medical Device Registered FDA Class 1 Medical Device Basic UDI-DI: 0710859306009
<p><u>SPECIFICATIONS:</u></p> <ul style="list-style-type: none"> OD 7 from 315 -11,000nm 315 nm - 1,400 nm RI LB 7, D LB 6 1,400 nm - 11,000 nm DI LB 4, R LB 3 Contact us for the latest test results from our independent Laboratory. 	<p><u>PACKAGING:</u></p> <ul style="list-style-type: none"> 25 pair per box Inner Box Dimensions: 2.75 in x 3.38 in x 1.62 in, 7 cm x 8.6 cm x 4.1 cm Master Box Dimensions 18inx12inx10in, 45.7 cm x 30.5 cm x 25.4cm 100 boxes per master carton, 25 pair per Box Master Carton - 14lbs, 6.35kg Box 25 pair - 2.2oz, 62gm



REF TT63

0 710859 306009

THEIA EYE BLOCK
3 Canal St.
Cumberland RI 02864
USA - Tel (888) 626 3319
www.theiaeye.com

Single Use

THEIA EYE BLOCK

LASER EYE SHIELD - 25PR/PK

REF TT63 **43877-1-1** **LOT**

THE-Laser-G2 Rev. C

0 710859 306009

QTY - 100 PK

100 PK

Jan/2023

0710859306061

Made in the U.S.A.