



Walmay Medical Products

presents the new and improved

CRYOSPRAY

**It is the instrument of choice for all
Dermatologists and Practitioners.**

- Ergonomic collar design
- 350ml Flask – easy for all Operators to hold.
- Fixed Stainless Steel Probe with greater tensile strength
- 5 tips, 4 Coned for accurate aim during procedure and 1 curved tip for body orifices.
- Regulator Handle may be placed on either side of the probe to suit left or right hand Operator
- Centre of gravity is towards the top of the flask, making it easy to balance and steady the instrument during use.
- Stainless Steel and Brass components for high quality instrument.
- Safety Valve maintains consistent pressure in the unit for operation.
- Carry Bag for safe and convenient transport.



**NEW &
IMPROVED
MODEL**

**5 YEAR
WARRANTY**

**MADE IN
AUSTRALIA**

docStock

MEDICAL EQUIPMENT SPECIALISTS™

DocStock Pty Ltd

14/11-11/15, 2015

P: (02) 8850 7177

E: info@docstock.com.au

www.docstock.com.au

W
Walmay Medical
PRODUCTS

Cryospray

WM227

For the controlled
application of
liquid nitrogen
in Cryosurgery.

ARTG No. 170073 C E 0805

Walmay Integrated Pty Ltd

P: +61 2 9432 2899

E: sales@cryospray.com.au

Made in Australia

REV 4