The SureTemp® Plus 690—

of care even higher

setting the standard



## SureTemp® Plus 690 specifications

....±0.2 °F (±0.1 °C) (Monitor Mode) Calibration Accuracy... Temperature Predict Time: ..Approx. 4-6 secs Adult Axillary (18 years and older)... ...Approx. 15 secs Pediatric Axillary (17 years and younger) ... ...Approx. 10 secs Rectal.... ...Approx. 10 secs Unit Dimensions ... ...8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm) ...12.6 oz (357 g) ...50.0 °F (10 °C) to 104 °F (40.0 °C) Ambient Temperature Range Patient Temperature Range... ..80.0 °F (26.7 °C) to 110 °F (43.3 °C) Power Source ...(3) 1.5 Vdc AA batteries Battery Life ..Approx. 6,000 temperatures Warranty ...2-year limited instrument, 1-year limited probe, and 90-day limited probe well

## ordering information

thermometers

01690-200 SureTemp Plus Electronic Thermometer (Model 690), Oral Probe with Oral Probe Well, 2-Year Limited Warranty

01690-201 SureTemp Plus Electronic Thermometer (Model 690), Rectal Probe with Rectal Probe Well, 2-Year Limited Warranty

01690-300 SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft Oral Probe with Oral Probe Well, 2-Year Limited Warranty (must be locked to the wall)

01690-301 SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft Rectal Probe with Rectal Probe Well, 2-Year Limited Warranty (must be locked to the wall)

01690-400 SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, Oral Probe with Oral Probe Well, 2-Year Limited Warranty

01690-700 SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, 9ft Oral Probe with Oral Probe Well, 2-Year Limited Warranty, with Rolling Stand and Holder for Spare Probe & Well accessories

21326-0000 Wall Holder

**21330-0000** Holder for Spare Probe and Well **06138-000** Calibration Key Assembly, 690/692

 02892-000
 Probe and Well Kit, 4 ft, Rectal

 02892-100
 Probe and Well Kit, 9 ft, Rectal

**02893-000** Probe and Well Kit, 4 ft, Oral **02893-100** Probe and Well Kit, 9 ft, Oral

**02891-0000** Probe Well, Blue **02891-1000** Probe Well, Red

08273-000 Thermometer Stand with Holder for Spare Probe and Well

21333-0000 Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well 98024-0000 Extended Warranty M690

replacement probe covers

**05031-101** Case of 1,000 Probe Covers

**WelchAllyn**°

Call Welch Allyn to order the SureTemp® Plus 690
4341 State Street Road, PO Box 220
Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 • (f) 315.685.3361
www.welchallyn.com
SureTemp is a registered trademark of Welch Allyn © 2006

Advancing Frontline Care™

SM2546 Rev D





## Introducing the SureTemp® Plus 690 Electronic Thermometer

Now taking temperatures is easier, safer, and more accurate. And taking care of all your patients is more convenient.



Advancing Frontline Care™



- Easy to use, intuitive operation
- Removable probe well reduces the risk of cross-contamination
- Improved accuracy for even more reliable readings
- Greater versatility for use with all your patients, all the time

The SureTemp® Plus 690—with the speed and accuracy you require, and the versatility you need for use with all your patients, all the time.

When we set out to improve the Welch Allyn SureTemp, we knew we had our work cut out for us. After all, when you start with the thermometer that helped set the standard for speed and accuracy, the expectations run pretty high.

So how did we do? See for yourself by taking a look at the innovative SureTemp Plus 690—easier, safer and more accurate than the original, and still built to withstand constant use in demanding environments.

Easy: An easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape make using the SureTemp Plus 690 a breeze.

**Safe:** Patient safety is improved thanks to an interchangeable probe well that reduces the risk of cross-contamination.

Fast and Accurate: In addition to improved accuracy, the SureTemp Plus 690 is one of the fastest thermometers on the market, providing approximately 4-6 second oral, 10-second pediatric axillary, 15-second adult axillary, and 10-second rectal measurements.

Versatile: To top it all off, you can now use the SureTemp Plus 690 with all your patients, from infants to adults—whether you need to take oral, axillary or rectal temperatures.



Interchangeable probe with



Battery life indicator

User-selectable icon and/or words for identification of thermometer mode

Mode button for selecting oral, pediatric axillary, or adult axillary





Interchangeable rectal probe and well reduce risk of cross-contamination



The SureTemp® Plus 690

