



Siemens Definition Flash 128

DATA SHEET

- **Brand:** Siemens
- **Model:** Definition Flash
- **YOM:** 2012
- **ID number:** SMS-R57984
- **Tube YOM:** 2022
- **Tube count:** 90.753 ss (as 27/04/2023)
- **Tube type:** Straton MX P
- **Tube installation date:** 13-10-2022
- **Tube 2 YOM:** 2020
- **Tube 2 count:** 206.637 ss (as 27/04/2023)
- **Tube 2 type:** Straton MX P
- **Tube 2 installation date:** 28-09-2020
- **Slices:** 128
- **Detector type:** Stella
- **Cooling:** Water cooling
- **SW version console:** Syngo CT VA48A
- **SW options console:** GET_WORKLIST, CAP3D_FILTER_SSD, CAP3D_FILTER_MIP, CAP3D_MAIN, CAP3D_VOIMODE, CARDIO_PRO, CARDIO_RETRO, BOLUS_TRACK, CARE_DOSE, SYNGO_GENERAL_LICENSE, CAP3D_FILTER_VRT, CAP3D_EDITOR, MPPS_SUPPORT, VOLUME_WORKFLOW, CTAPP_MAIN, DUAL_MON, INJ_CONNECTION, UHR_03, UHR, CAP3D_BONEREMOVAL, CAP3D_MULTIMONITOR, REMOTE_DESKTOP, DICOM_SR_VIEWER, CTAPP_FILTER_VRT, CTAPP_FILTER_MIP, CTAPP_FUSIONBASE, CTAPP_FUSIONCLIPVIEW, CTAPP_FUSIONLAYOUT, CARDIO_BESTPHASE_FINDING,



DUAL_ENERGY_SCAN_RECON, CTAPP_FILTER_SSD,
FAST_3D_ALIGN, CT_NEURO_DSA, CTAPP_VPCT_NEURO,
CTAPP_CVF_PHASEVIEWER, VOL_PERFUSION_SCAN,
MAX_POWER, INJ_CONNECTION_ROW, SCAN_AUTOREGION,
CARE_DOSE_ADJUST, CARE_DOSE_PROFILE,
CARE_DOSE_DASHBOARD, RECON_ITERATIVE_STEP2,
CTAPP_DUALENERGY_MONO_E, DEFINITION_FLASH,
CARDIO_HALFBEAT, SPIRAL_HIGHPITCH,
CARE_DOSE_X_CARE, VOLTAGE_70KV,
CARE_DOSE_AUTOKV, SCAN_GUIDE, RECON_SPINEREGION,
RECON_ITERATIVE_STEP1, DUAL_ENERGY_FAST,
0.5_MM_SLICE, MR_SUPPORT_i,





