


PRO-E Series of EEG Needle Electrodes

When the first few minutes can make a life-long difference



Reduce the risk of needle sticks – grasp the cap and pull the lead-wire to retract the needle!

Key Features

- When time is of the essence - for neonatal/NICU, pediatric, or adult diagnostic EEG/aEEG
- Unique Protectrode retractable blue cap protects you from needle sticks
- 12mm, 29 gauge (0.29 mm diameter) uncoated medical grade stainless steel needles
- Good signal quality with low impedance
- Very low-profile hub allows the needle to lay flat against the scalp or skin
- Conical tip is designed to be sharp and easy to insert with no hooked tips and minimal bleeding
- Available in assorted lead-wire colors, lead-wire lengths and package quantities
- Standard 1.5mm touch-proof DIN 42802-1 safety connectors
- Sterile; 3-year shelf life



PRO-E4 or PRO-E5 (2 purple, 2 black); PRO-E3 (red, yellow, black); PRO-E12SAF (10 color set)