


Disposable Monopolar EMG Needle Electrode - pencil point



Technomed understands the importance of patient comfort and signal quality in EMG.

Technomed's Monopolar Needle Electrode - pencil point and cable solutions are designed for optimal balance between needle sharpness and control of movement inside the muscle.


The hub and cable provide a firm grip and optimal precision in needle positioning.

User Benefits

- Clear and reliable signal
- Optimal patient comfort with ease of penetration and repositioning
- Tip design reduces bleeding
- Precise handling with improved ergonomic hub and cable design
- Approved for stimulation
- Optional shielded cable for 50/60 Hz interference elimination

Available variants

Product Description	Needle Length	Diameter	Hub Color	Part Number
Pre-wired	25 mm	0.35 mm (28 Gauge)	Red	TE/AP2535-335
Pre-wired	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AP2545-335
Pre-wired	37 mm	0.35 mm (28 Gauge)	Orange	TE/AP3735-335
Pre-wired	37 mm	0.45 mm (26 Gauge)	Green	TE/AP3745-335
Pre-wired	50 mm	0.45 mm (26 Gauge)	Blue	TE/AP5045-335
Pre-wired	75 mm	0.45 mm (26 Gauge)	Violet	TE/AP7545-335
Detachable	25 mm	0.35 mm (28 Gauge)	Red	TE/AD2535-008
Detachable	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AD2545-008
Detachable	37 mm	0.35 mm (28 Gauge)	Orange	TE/AD3735-008
Detachable	37 mm	0.45 mm (26 Gauge)	Green	TE/AD3745-008
Detachable	50 mm	0.45 mm (26 Gauge)	Blue	TE/AD5045-008
Detachable	75 mm	0.45 mm (26 Gauge)	Violet	TE/AD7545-008



Packing: 25 individually packed electrodes per box
Pre-wired: 0.75 m

Reusable Cable	Wire Length	Part Number
Shielded	1 m	TE/D54-451
Unshielded	0.6 m	TE/D34-231

Specifications

Lead-free	✓
Latex-free	✓
BPA-free	✓
DEHP < 0.1% by weight	✓
REACH & RoHS compliant	✓
EtO sterilized	✓
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802
Shelf life for unopened pouches	36 months



Your Technomed Business Partner:



The Electrode Store www.electrodestore.com
 360-829-0400 info@electrodestore.com



Taking you further. Step by step.

Technomed
 Amerikalaan 71
 6199 AE | Maastricht-Airport | The Netherlands
 T +31 43 408 68 68 | info@technomed.nl

Technomed USA
 Suite 4990 Constellation Drive
 White Bear Township | MN 55127 | United States
 T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl