



Good-Lite Company
1155 Jansen Farm Drive
Elgin, IL 60123

good-lite.com
800.362.3860
orders@good-lite.com

Surface cleaning of diagnostic eye products (including Vision Tests such as Wall charts and Near Vision cards, Stereo Tests, Trial Frames, Frosted Lens Occluders, Confirmation flippers, Red/Green Glasses, Prisms, Trial Lenses)

Cleaning Guide

- Use non-solvent-based germicide, such as alcohol, non- alcohol, or antiseptic disinfecting wipes.
- The product can be carefully wiped with a spray-dampened lint-free cloth.
- Do not use harsh solvents, such as acetone.
- Do not use microfiber cloths.