

Isolation Gown Level 1, DG-5

Isolation Gown - Level I	
AAMI Level / Performance Requirements At 4% AQL	Spray Impact Penetration <= 4.5g
Anticipated Risk of Exposure to Fluid	Low
Gown Characteristics	Lightweight PP, Non-Woven
Recommended Areas of Use	Laundry, Services, Housekeeping



In today's environment, it is more important than ever to be equipped with the proper protective apparel to help reduce the spread of infections. AAMI guidelines are a widely accepted system of classification for protective apparel based on liquid barrier performance, and our Disposable Isolation Gown conforms to AAMI Level 1. Vizocare's Disposable Isolation Gown provides barrier protection based on AAMI guidelines. Our Isolation Gowns are made from PP Non-woven fabric and gives personnel the confidence and comfort they need to deliver.



Contact us:

32 W. 8th St., Suite No. 607, Erie, Pennsylvania 16501
 FDA Registration: 3017631823

Phone: (814) 840-1880
 Email: sales@gdimed.com

Isolation Gown Features:

- Breathable, flexible, and strong enough for tough duties
- Made from lightweight fluid-resistant multi-ply material
- Gowns are generously sized for full coverage and flexibility
- Offered with elastic or knit cuffs
- Constructed with waist and neck ties to ensure secure fit



Product Description:

Material: 25 GSM PP, Non-Woven

Type: High quality disposable nonwoven isolation dressing gowns

Certification: FDA, CE

Color: Blue, Yellow

Package Size: 16.5 x 10.6 x 12.2" In

Packaging: 10 Pieces per Bag, 100 Pieces per Case



Technical Specs:

Item: Disposable Isolation Gown

Protection Level: AAMI Level 1

Fabric: Lightweight Non-Woven

Size: Universal size, 47" X 55"

Style: Tie on Collar and Waist, Knit Cuff, Closed Front

Recommended usage: Basic care, Laundry, Housekeeping

