



Vizocare 3-Ply Face Mask for Kids, FM-7

Technical Specifications:

Type	Disposable face mask for kids
SKU	FM-7
Material	30% Melt-Blown fabric 70% Non-Woven fabric
Size	5.7 x 3.5 Inches
Color	Blue
Noise clip	Yes
Protection Class	Level 2
Fluid Resistance	120 mmHg
BFE	>98%
Certification	FDA, ISO
Applied Standard	GB15979-2002
Packaging	50 per box, 40 boxes per case, 30 cases per pallet

Features:

- Multi-Layer filtering effectively blocks dust and pollen in the air.
- Upgraded elastic ear band, ultra-comfortable to wear.
- Ultrasonic spot welding
- Elastic loops fit comfortably around ears
- Adjustable nose clip for excellent fit
- Fluid resistant and disposable
- Non-Sterile



Contact Us

Phone: (814) 840-1880
Email: sales@gdimed.com
www.gdimed.com

32 W. 8th St., Suite No.
607, Erie, Pennsylvania
16501
FDA Registration:
3017631823



(814) 840-1880
sales@gdimed.com
www.gdimed.com



32 W. 8th St., Suite No. 607
Erie, Pennsylvania 16501
FDA Registration: 3017631823

